

The port reservoirs without compromise

# NuPort®-CT

Implantable, venous, CECT compatible Port Catheter Systems



PHS Medical GmbH Ederweg 3 34277 Fuldabrück . Germany Tel. +49 561 998597-0 . Fax -199 info@phs-medical.de www.phs-medical.com

## NuPort®-CT

Implantable, venous, CECT compatible Port Catheter Systems

The most effective solution for long-term vascular access from specialists. In order to achieve optimal flushing and flow characteristics we have intensively studied fluid mechanics in the development of the Nu-Flow chamber. The result is an optimized port chamber geometry, which creates turbine-like flow dynamics. This compact, cosmetically inconspicuous port is especially suitable for cachectic or small patients and while incorporating all the advantages of the NuPort® family. The transparent septum also makes it possible to view the chamber during implantation in order to verify the venting




ISO 13485 Certified Quality System CE 0483

#### The benefits at a glance

- suitable for power injection (CECT compatible)
- · greatly reduced risk of occlusion and infection
- significantly improved flushing efficiency and minimal maintenance required
- facilitates blunt preparation of the port pocket
- MRI compatible up to 3 Tesla
- · a variety of catheter types

#### Plus simplified

- blood product administration & blood collection
- · parenteral nutrition
- · high viscosity medications



CECT and MRI compatible.

#### Port chamber

and flushing of the port. The tissue-compatible titanium is

The elliptical chamber, in combination with the tangential outlet, improves flow dynamics and reduces the risk of premature occlusion.



#### Detection

The NuPort®-CT is characterized by its special heart shape and the hexagonal septum opening. This is also clearly shown in the X-ray image. Furthermore, the heart shape supports the blunt preparation of the port pocket during implantation.



#### Connection

The connection mechanism enables, via audible and tactile feedback, a secure connection of the catheter to the port base. These two detection functions ensure a safe connection is guaranteed.



#### Localization

The raised edge of the NuPort®-CT allows easy palpation of the port and the the septum for infusion or CT application.

## Technical Data:

#### Material, volume, weight

| Port   | Titanium | Volume | ~ 0.5 ml |
|--------|----------|--------|----------|
| Septum | Silicone | Weight | 12.2 g   |

All catheters with rounded tip and length markings in 1 cm intervals up to 30 cm, radiopaque.

#### Flow rate and pressure in CECT application

| Needle size      | 19 Ga  | 20 Ga            | 22 Ga  |  |
|------------------|--------|------------------|--------|--|
| Flow rate        | 5 ml/s | 5 ml/s           | 2 ml/s |  |
| Maximum pressure | 300    | 300 psi/20.6 bar |        |  |


The max. flow rate or the pressure setting specified above must not be exceeded during power injection via the NuPort®-CT.

#### NuPort®-CT Basic System

Single lumen, venous system for implantation via "cut-down" technique.

#### Content:

NuPort®-CT, Catheter, Catheter lock, Huber needle, Blunt tip needle (flushing), Vein lifter (pick), 10 ml syringe.




## NuPort®-CT Complete Set

#### with percutaneous introducer set

Single lumen, venous system for implantation via percutaneous technique.

#### Content:

NuPort®-CT, Basic system, Introducer needle Split-Sheath cannula with dilator, Guidewire with "Thumb-Feed" guide, Tunneler



#### **Dimensions**





Ref. no. Basic

System w/o

**Introducer Set** 

CTP-066CP

CTP-075CP

CTP-008CS

CTP-009CP

CTP-096CS

### **Ordering Information**

| Cathatan           | C:    | 0.0     | I.D.    | Lamath |
|--------------------|-------|---------|---------|--------|
| Catheter           | Size  | 0.D.    | I.D.    | Length |
| Polyurethane       | 6.6 F | 2.20 mm | 1.20 mm | 70 cm  |
| Polyurethane       | 7.5 F | 2.50 mm | 1.15 mm | 50 cm  |
| Silastic® silicone | 8.0 F | 2.67 mm | 1.40 mm | 50 cm  |
| Polyurethane       | 9.0 F | 3.00 mm | 1.60 mm | 50 cm  |
| Silastic® silicone | 9.6 F | 3.18 mm | 1.58 mm | 50 cm  |

Ref. no. Complete
Set incl.
Introducer Set
CTP-066IP
CTP-075IP
CTP-008IS
CTP-009IP
CTP-096IS