PHS Medical GmbH

Ederweg 3
D 34277 Fuldabrück, Germany
Tel.: +49(0)561 998597-0
Mail: info@phs-medical.de

16.09.2024

Portsysteme inkl. Zubehör SSCP Rev. 02

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Busch, Jasmin D., et. al., Complication Rates Observed in Silicone and Polyurethane Catheters of TIVADs Implanted in the Upper Arm, *J Vasc Interv Radiol* 2017, 1–7

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Jordan, Karin, et. A., Venous Access Ports: Frequency and Management of Complications in Oncology Patients, *Onkologie* 2008;31:404–410

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Kunz-Virk, Julia, Krüger, Karsten, Power-injectable totally implantable venous access devices - analysis of success and complication rates of ultrasound-guided implantation, *Vasa* (2019), 48, 524–530

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Summary of Safety and Clinical Performance: NuPort®/C-Port® Port Systems Revision 02

Table of contents

Part 1: Information for users/health care professionals 3

This summary of safety and clinical performance summarizes the most important aspects regarding the safety and clinical performance of PHS Medical GmbH port systems. This summary does not replace the instructions for use to ensure safe use of the products.

Part 2: Information for patients/public 21

This Summary Safety and Clinical Performance Report (SSCP) is intended to provide the public with an updated summary of the most important aspects of the safety and clinical performance of the port systems manufactured by PHS Medical GmbH. The information below is intended for patients or lay persons. A more detailed summary of the safety and clinical performance of the product for healthcare professionals can be found in the first part of this document. The SSCP is not intended to provide general advice for the treatment of any medical condition. Please consult your healthcare professional if you have any questions about your medical condition or the use of the products in your situation. This SSCP is not intended as a substitute for the patient identification card or instructions for use, which contain information on the safe use of the products.

Appendix I: Library:..... 34

**Summary of Safety and Clinical Performance:
NuPort®/C-Port® Port Systems Revision 02**

Abbreviations and acronyms

Basic UDI-DI	Basic Unique device identification device identifier
BaSO ₄	Barium sulfate
CAPA	Corrective And Preventive Actions
CECT	Contrast Enhanced Computed Tomography
CER	Clinical Evaluation Report
CMR	Cancerogen Mutagen Reprotoxic, fertility-endangering
CVC	Central Venous Catheters
EtO	Ethylene oxide
F	French (specification of the outer diameter of catheters)
FSCA	Field Safety Corrective Action
G	Gauge
MDR	Medical Device Regulation
mPa	Milli Pascal
MRI (MRT)	Magnetic Resonance Imaging (<i>Magnetic Resonance Tomography</i>)
PICC	peripherally inserted central venous catheter
PMCF	Post Marked Clinical Follow-up
PMS	Post Marked Surveillance
PSI	Pounds per square inch (unit of measurement for pressure)
PUR	Polyurethane
RKL	Supplier complaint
SAR	Specific Absorption Rate
SSCP	Summary of Safety and Clinical Performance
U	Units

**Summary of Safety and Clinical Performance:
NuPort®/C-Port® Port Systems Revision 02**

Part 1: Information for users/health care professionals

1. Labeling of the product and general information

1.1. Trade name(s) of the device

NuPort®/C-Port®

Table 1: Article list titanium ports

Model	Description
CTP-XXXCP	NuPort-CT, port system set, standard, CT compatible, power injection capable, with PUR catheter
CTP-XXXIP	NuPort-CT, port system incl. introducer set, standard, CT compatible, power injection capable, with PUR catheter
CTP-XXXCS	NuPort-CT, port system set, standard, CT compatible, power injection capable, with silicone catheter
CTP-XXXIS	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
CTP-L-XXXCP	NuPort-CT, port system set, CT compatible, power injection capable, with PUR catheter, left chest placement
CTP-L-XXXIP	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with PUR catheter, left chest placement
CTP-L-XXXCS	NuPort-CT, port system set, CT compatible, power injection capable, with silicone catheter, left chest placement
CTP-L-XXXIS	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter, left chest placement
LPA-XXXCP	NuPort-LP, low-profile port system set, CT compatible, power injection capable, with PUR catheter
LPA-XXXIP	NuPort-LP, low-profile port system incl. introducer set, CT compatible, power injection capable, with PUR catheter
LPA-XXXCS	NuPort-LP, low-profile port system set, CT compatible, power injection capable, with silicone catheter
LPA-XXXIS	NuPort-LP, low-profile port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
SPB-XXXCP	C-Port, port system set, standard, with PUR catheter
SPB-XXXIP	C-Port, port system incl. introducer set, standard, with PUR catheter
SPB-XXXCS	C-Port, port system set, standard, with silicone catheter
SPB-XXXIS	C-Port, port system incl. introducer set, standard, with silicone catheter
PH-XXIS-YY	NuPort-CT, port system incl. introducer set, standard, with silicone catheter
LPA-XXXIP-YY	NuPort-LPA, low-profile port system incl. introducer set, with PUR catheter
CTP-XXXCP-YYY	NuPort-CT, port system set, CT compatible, with PUR catheter
CTP-XXXIP-YYY	NuPort-CT, port system incl. introducer set, CT compatible, with PUR catheter
CTP-L-XXXCP-YYY	NuPort-CT, port system set, CT compatible, with PUR catheter, left chest placement

Summary of Safety and Clinical Performance: NuPort®/C-Port® Port Systems Revision 02

CTP-L-XXXIP-YYY	NuPort-CT, port system incl. introducer set, CT compatible, with PUR catheter, left chest placement
8086551	Fresenius C-Port HF 8F, 8F silicone catheter, with introducer set
8086561	Fresenius C-Port HF 9.6F, 9.6F silicone catheter, with introducer set

Table 2: Item list PSU ports

Model	Description
CTKP-XXXCP	C-Port-CT, port system set, standard, CT compatible, power injection capable, with PUR catheter
CTKP-XXXIP	C-Port-CT, port system incl. introducer set, standard, CT compatible, power injection capable, with PUR catheter
CTKP-XXXCS	C-Port-CT, port system set, standard, CT compatible, power injection capable, with silicone catheter
CTKP-XXXIS	C-Port-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
KP-XXXCP	C-Port-II, port system set, standard, with PUR catheter
KP-XXXIP	C-Port-II, port system incl. introducer set, standard, with PUR catheter
KP-XXXCS	C-Port-II, port system set, standard, with silicone catheter
KP-XXXIS	C-Port-II, port system incl. introducer set, standard, with silicone catheter
CTKP-XXXIS-YY	C-Port-CT, port system incl. introducer set, standard, with silicone catheter
CTKP-XXXIP-YY	C-Port-CT, port system incl. introducer set, standard, with PUR catheter
CTKP-XXXCP-YYY	C-Port-CT, port system set, CT compatible, with PUR catheter

XXX - catheter size

YY/Y - abbreviation for customers of private label port systems (distribution by customers, placing on the market by PHS Medical GmbH)

The ports are available with catheters in different sizes. The following catheters are available:

- 5F PUR
- 6.6F PUR
- 7.5F PUR
- 9F PUR
- 8F Silicone
- 9.6F Silicone

1.2. Name and address of the manufacturer

PHS Medical GmbH
Ederweg 3
34277 Fuldabrück
Germany

Tel.: +49 (0) 561 - 99 85 97 0

E-mail: info@phs-medical.de

**Summary of Safety and Clinical Performance:
NuPort®/C-Port® Port Systems Revision 02**

1.3. Single registration number (SRN) of the manufacturer

DE-MF-000005179

1.4. Basic UDI-DI

Titan ports: 42503746PortsystemTitanNF

PSU ports: 42503746PortsystemPSUC6

1.5. Nomenclature for medical devices

"Port Systems" 61494

1.6. Product class

Class III

(Regulation (EU) 2017/745 Annex VIII Rule 8)

*"All implantable devices and long-term surgically invasive devices are classified as class IIb unless...
... they are intended to be used in direct contact with the heart, the central circulatory system or the central
nervous system, in which case they are classified as class III."*

1.7. Year in which the first certificate (CE) for the product was issued

2003

1.8. Authorized representative, if applicable; name and SRN

Not applicable.

1.9. Notified body

mdc medical device certification GmbH
Kriegerstrasse 6
70191 Stuttgart
Germany

Tel.: +49 (0) 253597-0

E-mail: mdc@mdc-ce.de

Identification number 0483

**Summary of Safety and Clinical Performance:
NuPort®/C-Port® Port Systems Revision 02**

2. Intended purpose of the product

2.1. Intended use

NuPort®/C-Port® port systems are subcutaneously implantable medical devices used when repeated access to the vascular system is indicated as part of therapy. These products provide the user with an easy-to-find site for cannula insertion to administer fluids, such as drugs, into the vascular system or to draw blood samples. The port is accessed by percutaneous puncture with a (non-coring) port cannula.

The **NuPort®/C-Port®** port systems are available as standard vascular access ports NOT suitable for power injections and as port systems for CECT high pressure injections. Power injections are performed using an appropriate infusion set with a (non-coring) port cannula. Port systems are intended for single use and are explanted and disposed of at the end of therapy.

2.2. Indications and target group

The **NuPort®/C-Port®** port systems are intended for the following patient uses:

- repeated access to the vascular system for injections
- drug infusion
- administration of blood or blood products and/or collection of blood as part of therapy.

The **NuPort®/C-Port®** port systems can be implanted in different areas of the body.

Areas of application include:

- Parenteral nutrition
- Chemotherapy
- Pain therapy
- endocrinological therapy
- Substitution of coagulation factors
- Contrast medium infusion in radiation diagnostics
- Therapy of severe asthma or cardiac arrhythmias
- Blood collection and transfusion

NuPort®/C-Port® are intended for all patients requiring permanent vascular access. Port and catheter size should be selected according to the patient's age and stature.

2.3. Contraindication and restrictions

The **NuPort®/C-Port®** port systems are contraindicated under the following conditions:

- Patients with known or suspected allergic reactions to the materials contained in the implant
- Patients with known infection, bacteremia, or sepsis
- Patients with body tissue that cannot adequately support the port or catheter
- Patients with chronic obstructive pulmonary disease (COPD)
- When a potential access site has been irradiated
- If the patient's tissue factors prevent proper stabilization of the product and/or access
- Previous venous thrombosis or surgical procedures at or near the site of the planned use

**Summary of Safety and Clinical Performance:
NuPort®/C-Port® Port Systems Revision 02**

3. Product description

3.1. Description of the medical device

The product is a set of various components which is sterilized using EtO (ethylene oxide). It is a single-use product and must not be reused under any circumstances.

The **implantable** parts include a port, which is made of titanium or polysulfone (PSU) depending on the variant, and in which a silicone septum is mounted, as well as a catheter made of silicone or polyurethane (PUR) and a catheter fixation, which is made of polysulfone or polycarbonate. The septum is designed to withstand multiple punctures with a non-coring port cannula. To connect the catheter to the port, a connector is screwed into the port and bonded with silicone adhesive. The catheter is slid onto this during implantation and fixed with the catheter fixation for a secure connection.

The lifetime of the implantable parts is limited by the number of punctures of the septum. Tests have shown that 1,000 punctures/cm² are possible when using a 22-gauge Huber cannula. After this number of punctures, the septum may leak. These values may vary individually depending on the size of the cannula used and the technique of drug delivery. In the interest of patient safety, the possibility of leakage and elective replacement should be considered after 1,000 punctures.

For implantation, the following components for **single use** are also included in the set:

- Huber cannula, which enables puncture of the septum without coring
(Material: cannula - stainless steel, luer - polycarbonate)
- Blunt cannula for flushing the catheter
(Material: cannula - stainless steel, luer - polycarbonate)
- 10ml syringe with luer slip
(Material: Polypropylene, Polyisoprene)
- Vein pick, which is intended to facilitate the insertion of the catheter into the vein
(Material: Terlux)

Furthermore, the portsets are also available as a variant with introducer. This is used for insertion of the catheters via the minimally invasive percutaneous technique (cf. "Seldinger technique") or in the so-called open surgical technique.

The following additional components are included in the set with Introducer:

- Puncture cannula (Seldinger cannula)
(Material: cannula - stainless steel, luer - polycarbonate)
- Guide wire with dispenser and insertion aid
(Material: stainless steel, polypropylene, polyethylene)
- Split sheath (tearaway) with dilator
(Material: polyflon, polypropylene, polyethylene)
- Tunneler
(Material: stainless steel)

**Summary of Safety and Clinical Performance:
NuPort®/C-Port® Port Systems Revision 02**


Fig 1: Portset with Titanium port and Introducer

The **NuPort®/C-Port®** port systems do **not** contain **any**

- Drug substances (including derivatives from human blood or plasma),
- Tissues or cells of human or animal origin or their derivatives,
- Substances or combinations of substances that are absorbed by the human body or distributed locally in the body,
- Materials containing CMR substances (carcinogens, mutagens or substances toxic to reproduction) or endocrine disruptors.

Although silicone is a known allergen, there is no evidence of allergies to medical silicone in the literature or from practice. There are no other known allergies or intolerances for the materials used.

3.2. Earlier variants/generations

There has been no change to the product since 2014. In 2014, the F6.6 PUR catheter was added.

3.3. Description of accessories

As an accessory, the Introducer Set (Basic-UDI: 42503746IntroducerYF) is available as a separate medical device, giving the operating physician the option of choosing between Seldinger or cut-down technique.

Contents of the separately available percutaneous introducer set

- Puncture cannula (Seldinger cannula)
- 10ml syringe with luer slip

**Summary of Safety and Clinical Performance:
NuPort®/C-Port® Port Systems Revision 02**

- Guide wire with dispenser and insertion aid
- Split sheath (tearaway) with dilator
- Tunneler

3.4. Description of other medical devices that can be used together with the product

The administration of drugs or parenteral nutrition solutions can be carried out by means of an infusion pump via the port systems. Infusion pumps are dosing pumps which are used for the continuous administration of infusions.

In the case of port systems suitable for CECT high-pressure injection, the contrast medium can be administered for radiation diagnostics by means of a high-pressure injection with a pump suitable for this purpose.

For the exact description of these pumps, please contact the pump manufacturer directly.

4. Risks and undesirable effects

4.1. Residual risks and undesirable effects

Use of the **NuPort®/C-Port®** family of ports involves risks commonly associated with the insertion or use of an implant or indwelling catheter, including but not limited to the following:

- Infection
- Erosion of the skin
- Hematoma
- Necrosis
- Displacement of the product
- Clot formation
- Catheter fragmentation
- Embolization
- Catheter closure
- Cardiac arrhythmias
- Cardiac puncture
- Pericardial tamponade
- Extravasation of drugs
- Fibrin sheath
- Endocarditis
- Implant rejections
- Port/Catheter Migration
- Perforation or rupture of vessels
- Injury of the thoracic duct
- Thromboembolism
- Thrombosis
- Superior vena cava syndrome
- Pneumothorax
- spontaneous malpositioning or retraction of the catheter tip
- Vascular erosion

Summary of Safety and Clinical Performance: NuPort®/C-Port® Port Systems Revision 02

- Catheter- or port-related sepsis
- Air embolism
- Bleeding
- Catheter embolism

The following table shows the most frequent complications reported to PHS Medical GmbH in the last 20 years and their frequency in scientific literature and from studies performed on port systems of PHS Medical GmbH:

Table 3: Frequency of complications of port systems

Complication	Average Complication rates from the literature*	Average complication rates from studies performed on products of PHS Medical GmbH	Number of complications reported to PHS Medical GmbH (2002-2022)	Percentage share in relation to the total volume of products sold to date at PHS Medical GmbH (2002- 2022)
Complications due to disregard of the instructions for use	Not named	0,5%	39	0,011%
"Pinch off", catheter rupture and dislocation	1%	0,9%	29	0,013%
Damaged septum	Not named	0%	27	0,012%
Occlusion	0,9%	0,9%	2	0,001%
Infection	4,5%	4,2%	1	0,0003%
Thromboses	2,8%	0,9%	0	0%

4.2. Warnings and precautions

Warnings:

- **Do not reuse!**
The products from PHS Medical GmbH must not be reused. Failure to observe this warning may cause allergic reactions or infections/cross-contamination in the patient.
- **Do not resterilize!**
The products from PHS Medical GmbH must not be re-sterilized. Failure to observe this warning may cause allergic reactions or infections/cross-contamination for the patient.
- **Do not use alcohol, acetone or solutions containing these substances**
These solutions can affect the properties of the plastic components, resulting in damage to the product.
- **No use by non-professionals**
Implantation and application may only be performed by trained personnel, otherwise serious handling errors may occur.
- **"Pinch off"**

**Summary of Safety and Clinical Performance:
NuPort®/C-Port® Port Systems Revision 02**

- The so-called "pinch off" occurs when the catheter is compressed between the clavicle and the first rib (costoclavicular space), especially when accessing via the subclavian vein.
- **No injection through a port with signs of pinch-off**
DO NOT attempt to inject through a port that shows signs of "pinch-off" as this may cause the port system to fail.
 - **Avoid damaging the catheter during implantation**
During implantation, care must be taken not to constrict the catheter during suturing or damage it with a cannula or scalpel, as the catheter can tear quickly once damaged.
 - **Do not use heparin in heparin-sensitive patients or patients with heparin-induced thrombocytopenia**
Some patients may be hypersensitive to heparin or have thrombocytopenia (HIT). The use of heparin in these patients is contraindicated.
 - **High-pressure injection only with products suitable for this purpose**
DO NOT attempt power injection through a port that is **not** labeled "CT compatible / Power Injectable".
 - **Bending of the cannula tip in case of excessive pressure**
The pressure exerted on the cannula by contact with the port floor leads in some cases to a bending ("fishhook") of the cannula tip. The size of the "fishhook" depends on the pressure applied by the user during puncture, but is also strongly influenced by the Huber cannula used (size and design). This can lead to inappropriate damage to the septum (especially with the NuPort®-LP-CT (LPA)) and extravasation. Therefore, excessive pressure on the cannula should be avoided at all costs.
 - **Not more than 300 PSI / ≈21 bar for power injection capable ports**
When performing power injection through the **NuPort®/C-Port®** port system, do not exceed the pressure limit of 300 psi / ≈21 bar or the maximum flow rate setting on the power injection injector.
 - **Cleaning blocked catheters**
If you encounter resistance during injection or infusion, stop the procedure immediately. The resistance may be due to a blocked catheter and excess pressure in a blocked catheter may cause the catheter or port to rupture.
Blood clotting with subsequent occlusion of the catheter or port is possible if the port is not properly flushed. If you find that the catheter or port is blocked, consider flushing with an appropriate lysing agent such as urokinase. While a dose of 5,000 U is recommended, please refer to the drug's packaging for instructions and possible complications associated with the use of Urokinase.
 - **Lifetime of the port (septum)**
The septum has a limited service life. Tests have been performed which show that 1000 punctures/cm² are possible when using a 22-gauge Huber cannula. After this number of punctures, the septum may leak. Your experience may vary depending on the size of cannula you use and your technique of drug delivery. In the interest of patient safety, after 1,000 punctures, look for the possibility of leakage and elective replacement.
 - **MRI related heating**
In non-clinical testing, the **NuPort®/C-Port®** port systems produced the following maximum temperature rise during MRI examinations performed with a 15-minute scan (i.e., per pulse sequence) in a 3 Tesla MR system (3 Tesla /128 MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI):

Whole-body averaged SAR specified for the MR system:	2.9 Watt/kg
Calorimetry readings, whole-body averaged SAR:	2.7 watts/kg
Highest temperature change:	≤ 2.4°C*

**Summary of Safety and Clinical Performance:
NuPort®/C-Port® Port Systems Revision 02**

Temperature scaled to a whole-body averaged SAR of 4 W/kg: $\leq 3.6^{\circ}\text{C}^*$.

Local SAR at the position of the reference probe: 3.2 W/kg.

(*Based on worst-case, C-port DL)

- **Guide wire not compatible with MRI examinations**

The guide wire as part of the insertion set is not compatible with MRI examinations. This can lead to overheating or movement of the guide wire.

- **Compatibility of the drugs used**

PHS Medical GmbH cannot guarantee the compatibility of the drugs used with the various port systems.

- **Avoid injecting drugs into the subcutaneous tissue**

The cannula must be held securely at the base of the port below the septum during these procedures to avoid injecting the drug into the subcutaneous tissue. Do not make any angular movement and do not twist the cannula and syringe once they are in the septum.

- **Do not use damaged guide wires**

DO NOT attempt to straighten a kinked or bent guide wire.

If the guide wire is deformed during the procedure, it must be discarded. Use a new guide wire for implantation.

- **No change of the product**

Do NOT modify the product.

- **No retraction of the guide wire via the puncture cannula**

NEVER pull the guide wire back over the positioned puncture cannula! This may damage the guide wire!

Precautions:

- Sale only to physicians or their authorized representatives
- For the low-profile port (NuPort®-LP-CT (LPA)), use **only 20 G and 22 G** port cannulas for standard applications; use **only 22 G port cannulas** for high-pressure injections. Using a too large cannula when puncturing the LPA port may damage the septum.
- Measures to prevent "pinch off":
 - Catheters placed over the subclavian vein should be inserted into the vessel at the junction between the lateral and middle thirds of the clavicle. If the catheter is inserted more medially, the catheter may be compressed between the first rib and the clavicle. This in turn could result in damage to the catheter.
 - Control position of the catheter after implantation by imaging to ensure that the catheter is not displaced between the first rib and the clavicle.
 - When using access via the internal jugular vein, pinch off is almost completely avoided.
- Syringes smaller than 10 ml should not be used. The use of syringes smaller than 10 ml results in a higher pressure at which it cannot be guaranteed that the port or catheter will not be damaged.
- The cannula must be held securely at the base of the port below the septum when puncturing the port to avoid injecting the drug into the subcutaneous tissue. Do not make any angular movement and do not twist the cannula and syringe once they are in the septum.
- After each drug injection, flush the port with 10 ml of saline solution. Do not aspirate or run fluids retrograde without flushing the port. Standing blood in the catheter may cause occlusion of the catheter. Flushing is essential for further use of the port.
- Ensure that the viscosity of the contrast medium is 11.8 mPa or less prior to injection (this may require heating the medium to achieve the correct viscosity).

**Summary of Safety and Clinical Performance:
NuPort®/C-Port® Port Systems Revision 02**

- Before implantation, visually inspect the port. Do not use if holes, cracks or surface contamination are detected.
- Purge all air from the port and catheter before inserting the product into the patient.
- Be careful when placing the NuPort®-LP-CT (LPA) as an arm port: Insert the catheter ONLY into the veins of the upper arm! The port can be implanted in the forearm and the catheter tunneled there. This implantation technique reduces the risk of crushing the catheter in the elbow when moving the arm.
- Position the reservoir pocket so that the suture line does not run directly over the port. The suture line can be lateral, medial, superior, or inferior to the port body, and the pocket should be large enough to prevent stretching or trauma to the suture line. The port should not be placed too deep so as not to complicate localization, and not too shallow so as to cause erosion in the form of pressure necrosis. A depth of about 5 mm below the skin surface is recommended as the optimal placement depth.
- PHS Medical GmbH recommends flushing the **NuPort®/C-Port®** port systems with saline solution every three months during longer breaks in use to prevent occlusion and to ensure that the port is still functional. However, the frequency of flushing is a decision based on the type of application. Here, flushing intervals of once a week, once a month, every three months, or even longer intervals are common in practice. Taking into account that each port puncture is also associated with the risk of germ introduction, the intervals should be determined individually.
- The **NuPort®/C-Port®** port systems have been tested for safety in an MRI environment at the University of Southern California Hospital MRI Center in accordance with the American Society for Testing and Materials (ASTM). Both ports of the **NuPort®/C-Port®** port systems were found to be MR capable.
 - The devices can be safely scanned immediately after insertion under the following conditions:
 - Static magnetic field of 3 Tesla or less
 - Maximum spatial gradient magnetic field of 12,000 Gauss/cm (extrapolated) or less.
 - Maximum Specific Absorption Rate (SAR) reported by the MR system averaged over the whole body of 4 W/kg for 15 minutes of scanning (i.e. per pulse train).
 - First stage-controlled operating mode for the MR system
- The quality of MR images may be compromised if the area of interest is in the same area or relatively close to the position of the port. Therefore, it may be necessary to optimize MR imaging parameters for the presence of the implant.
- Fluoroscopy is recommended to verify proper placement of the catheter tip in the superior vena cava (ideally between the 5th and 6th rib).
- Attach the port to the underlying fascia with at least three independent sutures arranged around the base. Pass each suture through the base of the reservoir and the fascia. This prevents migration of the port and eliminates the possibility of it tipping over in the port pocket.
- To avoid excessive damage to the port by punctures with Huber cannulas, we recommend:
 - avoid excessive pressure on the cannula during access
 - not puncturing the port at the same septal site each time
- During venipuncture, place the patient in an inclined position, head down, to prevent the occurrence of air embolisms. (Trendelenburg positioning)
- The guide wire should be routinely inspected before use and discarded if the guide wire is deformed.
- If resistance is encountered when advancing or retracting the guide wire, determine the cause by fluoroscopy and correct it before proceeding with the procedure. DO NOT rotate the guide wire if

**Summary of Safety and Clinical Performance:
NuPort®/C-Port® Port Systems Revision 02**

significant resistance is felt. Since guide wires are delicate and fragile, special care must be taken when handling them.

- Check the locking of the dilator with the split sheath (45° bayonet lock) before use!
- After removing the dilator from the split sheath, immediately close the opening with a fingertip to prevent blood loss or air aspiration.
- Avoid damaging the skin, bones and surrounding organs when using the tunneler!
- Avoid kinks, bruises and injuries within the later used catheter section.

Note

The user is obliged to report all serious incidents related to the product to the manufacturer and the competent authorities in the country where the user is located.

4.3. Other relevant safety aspects, including a summary of all field safety corrective actions (FSCA including FSN), if applicable.

In 2017, a recall occurred due to the lack of X-ray visibility of the F9.6 silicone catheters. This incident was registered at the Federal Institute for Drugs and Medical Devices under case number 01229/17, and at PHS Medical GmbH under the internal complaint number RKL036 and in CAPA002. The processing of the incident was completed on 07 April 2017. All products containing affected catheters have been recalled.

The cause of the event was insufficient mixing of the materials. Although the base material (Q7-4765) and BaSO₄ were correctly weighed and measured, the mixture - according to the operator - was not thoroughly processed with the intended 25 cycles in the grinder. Therefore, the BaSO₄ distribution was not homogeneous.

A CAPA was opened at the catheter supplier and the relevant work instructions were adapted to ensure that the materials are sufficiently mixed together.

In the meantime, the supplier for the silicone catheters has been changed.

5. Summary of post-market clinical follow-up (PMCF)

5.1. Summary of clinical data on an equivalent product, if applicable

Not applicable. As a separate clinical evaluation has been prepared for the **NuPort®/C-Port®** port systems, which is based on own clinical data.

5.2. Summary of clinical data from testing of the device conducted prior to CE marking, if applicable

Not applicable. Since the **NuPort®/C-Port®** port systems have already been on the market for more than 20 years.

5.3. Summary of clinical data from other sources, if applicable

Port systems and introducer sets from PHS Medical GmbH were used in clinical studies at the Universities of Münster and Essen. In the study at the University of Essen by Schenck et al. (2012), 347 patients were

**Summary of Safety and Clinical Performance:
NuPort®/C-Port® Port Systems Revision 02**

implanted with **NuPort®** (titanium/silicone) from PHS Medical GmbH between December 1999 and June 2011. In most cases (91.7%), implantation was due to planned chemotherapy. The other applications were artificial nutrition (4.6%) and lack of suitable peripheral vascular access (1.4%) or pain therapy (2.3%). The implantations of the port systems could be performed without complications in all cases. Early explantation of the port system within the first 30 days was necessary in six cases, three of which were due to catheter dislocation and one each to port-associated septicemia, port pocket infection, and arm vein thrombosis. The average duration of stay was 239 days (2 - 2604 days), demonstrating that the port systems from PHS Medical GmbH are suitable to provide permanent (> 30 days according to MDR) access to the vascular system.

This is also confirmed by the study of the University of Münster by Becker et al. (2021), who implanted a port system from PHS Medical GmbH in 500 patients and followed the patients over a period of 12 months. Implantation by open Seldinger technique using the introducer set from PHS Medical GmbH was mostly free of complications (4 complications 0.8%) and only in one case the goal of providing permanent (> 30 days) vascular access was not achieved, as here the port system had to be explanted prematurely due to infection. The complications that occurred in these two studies correspond in type, number and severity to the adverse events of port systems already described in the literature and listed in the reporting registers. However, the actual incidents reported to PHS Medical GmbH are significantly lower than the numbers from these two studies, but this is not due to incomplete traceability of PHS Medical GmbH products, but rather to an inadequate reporting system or reporting awareness. A corresponding note has therefore been included in the instructions for use stating that serious incidents must be reported to both PHS Medical GmbH and the relevant authorities.

5.4. Overall summary of clinical performance and safety

Port systems and introducer sets have been used in humans for many years, are not subject to regular design changes, and the complications and side effects associated with their use are well documented. Introducer sets can be used to quickly and reliably insert appropriately sized catheters into the human venous and central venous circulatory system.

Central venous port catheter systems are becoming increasingly important nowadays, especially in oncology, as they provide permanent vascular access, significantly facilitate the administration of drugs, chemotherapeutics, high-caloric nutrient solutions and blood products (etc.) as well as the collection of blood, and minimize the risk of infection due to their subcutaneous location. One parameter for demonstrating the performance of the products from PHS Medical GmbH is therefore their ability to provide such permanent vascular access (according to MDR "Long-term" > 30 days). Here, a permanently functioning seal between the port chamber and the tissue (septum tightness) and the permanent transfer of fluids between the port chamber and the outlet site of the catheter (port-catheter connection and catheter tightness) should prevent irritation/damage of surrounding tissue and veins by the administered drugs - compared to irritation/damage when using other vascular access systems. The lifetime of port systems is limited by the number of punctures with one cannula (1000/cm²).

The reported complication rates from the literature are used as a basis for evaluating the clinical safety of the products from PHS Medical GmbH. Depending on the source, these range between 2.6 and 46 %. Taking into account the average number of complaints about the port systems and introducer sets from PHS Medical GmbH of less than 1 %, an overall internal complication rate of max. 2 % for the port systems incl. accessories and introducer sets was therefore determined or defined as acceptable.

Summary of Safety and Clinical Performance: NuPort®/C-Port® Port Systems Revision 02

Both the clinical studies on the port systems and introducer sets from PHS Medical GmbH (see 5.3. *Summary of clinical data from other sources*) and the overall consideration of the port systems and introducer sets from PHS Medical GmbH in the clinical evaluation show that the **NuPort®/C-Port®** port systems from PHS Medical GmbH under consideration provide safe and long-term access to a patient's vascular system for a variety of indications. The established parameters for demonstrating safety and performance are met by the products from PHS Medical GmbH.

As the design and manufacture of the products are compliant with the standards, the port systems and introducer sets currently meet the relevant General Safety and Performance Requirements of the MDR. According to the current state of knowledge, based on the state of the art as well as the product-specific results from tests, clinical data and scientific literature, the benefits outweigh the risks of using the port systems and introducer sets. The analysis and evaluation of potential risks has shown that there are no increased residual risks for patients, users or third parties within the scope of the intended use of the port systems and introducer sets, which can be confirmed by the product-related clinical data. The risk reduction measures are also appropriate.

The instructions for use appear appropriate with regard to the intended user group. It should be assumed that the products under consideration here will only be used and applied (healthcare and nursing staff) by healthcare professionals (surgeons or radiologists for port implantation or physicians for the introducer set). These users should be properly trained and informed about the complications associated with the product. These correspond to those commonly associated with the insertion or use of an implanted device or indwelling catheter and are appropriately known to the users. The residual risks associated with the use of the devices are acceptable in terms of the benefit to the patient.

The benefits of the **NuPort®/C-Port®** port systems and introducer sets for patients have been clearly demonstrated in the Clinical Evaluation Report (CER). Regular updates of the CER are part of the post-marketing surveillance and vigilance activities. If there are significant changes that impact the original data, the CER must be amended accordingly.

5.5. Ongoing or planned post-market clinical follow-up

PMS (Post Market Surveillance) plan and PMCF (Post Market Clinical Follow-up) plan and report for the **NuPort®/C-Port®** port systems are regularly prepared at least once a year. The PMS activities performed at PHS Medical GmbH include:

- Customer surveys on satisfaction as well as product and safety-related topics
- Review and evaluation of complaints and CAPAs
- Review and evaluation of observations, reports, events and incidents concerning products of PHS Medical GmbH
- Audit of risk management
- Inspection and evaluation of data from the production process
- Checking the instructions for use
- Review of relevant new and revised standards
- Review and evaluation of the relevant scientific literature
- Checking relevant databases and registers
- Recording of all data in the PHS Medical GmbH trend analysis, regular review of this data and the statistics contained therein.
- Review of the PMCF plan and PMCF report.

**Summary of Safety and Clinical Performance:
NuPort®/C-Port® Port Systems Revision 02**

The results are recorded and evaluated in the annual safety report and measures are defined if action is required.

This systematic, proactive market surveillance ensures that the quality, safety and performance of products are continuously reviewed throughout the entire product life cycle. This ensures that risks in the post-marketing phase are identified and remedied as early as possible, and that products are continuously improved so as not to endanger patients and users and to provide the best possible clinical benefit.

The evaluation performed in the safety report and PMCF Report has confirmed the safety and performance of the PHS port systems during their expected lifetime. No previously unknown adverse reactions, side effects, contraindications or emergent risks were identified during the PMCF period. It could be confirmed that the risk-benefit ratio is acceptable.

No possible systematic misuse or off-label use of the product was found, so that it could be shown herewith that the intended purpose of the product is correct.

6. Therapeutic alternatives

Only central venous catheters (CVC) and peripherally inserted central venous catheters (PICC) represent an alternative to port systems. However, port systems have the advantage over central venous catheters that they are less susceptible to infection, require less care, can be used for longer periods of time, have higher cosmetic acceptability, and are perceived as less burdensome by the patient. To prevent infection with a CVC, the insertion site must be cleaned daily. In addition, complications such as catheter rupture or thrombosis can also occur with the CVC.

7. User qualification and training

The **NuPort®/C-Port®** port systems may only be used by a physician or by medical personnel. Prior to each initial use/implantation, PHS Medical GmbH will provide product training.

8. Harmonized standards and applied common specifications

The **NuPort®/C-Port®** port systems and their accessories are manufactured in accordance with the applicable laws, standards, guidelines and regulations governing the manufacture of medical devices. Common specifications for our products are not available.

Table 3: List of standards

Standard	Designation
DIN EN 556-1: 2002-03 (Corrigendum 1:2006-12)	Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 1: Requirements for terminally sterilized medical devices
DIN EN 868-5: 2019-03	Packaging for terminally sterilized medical devices - Part 5: Sealable pouches and reels of porous materials and plastic film construction - Requirements and test methods
DIN EN ISO 5832-3: 2022-03	Implants for Surgery - Metallic Materials - Part 3: Wrought Titanium 6-Aluminum 4-Vanadium Alloy (Ti6-Al4-V)

**Summary of Safety and Clinical Performance:
NuPort®/C-Port® Port Systems Revision 02**

DIN EN ISO 7153-1: 2017-02	Surgical instruments - Metallic materials - Part 1: Stainless steel
DIN EN ISO 7864: 2016-12	Sterile hypodermic needles for single use
DIN EN ISO 7886-1: 2020-11	Sterile hypodermic syringes for single use- Part 1: Syringes for manual use
DIN EN ISO 10555-1: 2018-04	Sterile, single-use intravascular catheters -- Part 1: General requirements
DIN EN ISO 10555-6: 2020-02	Intravascular catheters - Sterile, single use catheters - Part 6: Subcutaneous implanted ports
DIN EN 10088-3: 2014-12	Stainless steels - Technical delivery conditions for semi-finished products, bars, wire rod, drawn wire, sections and bright products of corrosion resisting steels for general purposes.
DIN EN ISO 10993-1: 2021-05	Biologische Beurteilung von Medizinprodukten- Teil 1
DIN EN ISO 10993-3: 2015-02	Biological evaluation of medical devices - Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity
DIN EN ISO 10993-4: 2017-12	Biological evaluation of medical devices - Part 4: Selection of tests for interactions with blood
DIN EN ISO 10993-5: 2009-10	Biological evaluation of medical devices - Part 5: Tests for cytotoxicity: in vitro-methods
DIN EN ISO 10993-6: 2017-09	Biological evaluation of medical devices - Part 6: Tests for local effects after implantation
DIN EN ISO 10993-7: 2022-09	Ethylene oxide sterilization residuals
ISO 10993-10: 2021-11	Biological evaluation of medical devices - Part 10: Tests for skin sensitization
DIN EN ISO 10993-11: 2018-09	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity
DIN EN ISO 10993-12: 2021-08	Biological evaluation of medical devices- Part 12: Sample preparation and reference materials
DIN EN ISO 10993-18: 2021-03	Biological evaluation of medical devices - Part 18: Chemical characterization of medical device materials within a risk management process
DIN EN ISO 10993-23: 2021-10	Biological evaluation of medical devices - Part 23: Tests for skin irritation
DIN EN ISO 11070: 2019-04	Sterile single-use intravascular introducers, dilators and guidewires
DIN EN ISO 11135: 2020-04	Sterilization of health care products -- Ethylene oxide -- Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
DIN EN ISO 11138-2: 2017-08	Sterilisation of health care products Part 2: Biological Indicators for EtO sterilisation processes
DIN EN ISO 11607-1: 2020-05	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
DIN EN ISO 11607-2: 2020-05	Packaging for terminally sterilized medical devices -- Part 2: Validation requirements for forming, sealing and assembly processes
DIN EN ISO 11737-1: 2021-10	Part 1: Determination of a population of microorganisms on products
DIN EN ISO 11737-2: 2020-07	Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
DIN 13273-7: 2020-12	Determination of the x-ray attenuation of catheters - requirements for testing
DIN EN ISO 13485: 2021-12	Medical devices -- Quality management systems

**Summary of Safety and Clinical Performance:
NuPort®/C-Port® Port Systems Revision 02**

DIN EN ISO 14155: 2020-12	Clinical investigation of medical devices for human subjects - Good clinical practice
DIN EN ISO 14630: 2013-03	Non-active surgical implants-General requirements
DIN EN ISO 14644-1: 2016-06	Cleanrooms and associated controlled environments - Part 1: Classification of air cleanliness
DIN EN ISO 14644-2: 2016-05	Cleanrooms and associated controlled environments - Part 2: Monitoring to provide evidence of cleanroom performance related to air cleanliness by particle concentration
DIN EN ISO 14644-3: 2020-08	Cleanrooms and associated controlled environments - Part 3: Test methods
DIN EN ISO 14971: 2022-02	Medical devices -- Application of risk management to medical devices
ISO 15223-1: 2022-02	Medical devices -- Symbols to be used with medical device labels, labeling and information to be supplied
DIN EN 17141: 2021-02	Cleanrooms and associated controlled environments - Biocontamination control
DIN EN 62366-1: 2021-08	Medical devices -- Application of usability engineering to medical devices
DIN EN ISO 80369-7: 2021-08	Small-bore connectors for liquids and gases in healthcare applications - Part 7: Connectors for intravascular or hypodermic applications
ISO/TR 20416: 2020-07	Medical devices - Post-market surveillance for manufacturers
ISO/TR 24971: 2020-06	Medical devices - Guidance on the application of ISO 14971
DIN EN ISO 20417: 2022-03	Medical devices - Information to be supplied by the manufacturer
ISTA-2A: 2011	Simulation test for individual packaged-products less than 150 lbs
ASTM F88/F88M: 2021	Standard Test Method for Seal Strength of Flexible Barrier Materials
ASTM A276/A276M: 2017	Standard Specification for Stainless Steel Bars and Shapes
ASTM F136-13 (refer also to DIN EN ISO 5832-3)	Standard Specification for Wrought Titanium-6Aluminum-4Vanadium ELI (Extra Low Interstitial) Alloy for Surgical Implant Applications (UNS R56401)
ASTM F640-20	Standard Test Methods for Determining Radiopacity for Medical Use
ASTM F1929-15	Standard Test Methods for Detecting Seal Leaks In Porous Medical Packaging By Dye Penetration
ASTM F1886F/F1886M: 2016	Standard Test Method for Determining Integrity of Seals for Flexible Packaging by Visual Inspection
ASTM F1980-21	Standard Guide for Accelerated Aging of Sterile Barrier Systems for Medical Devices
ASTM F2052-15	Standard Test Method for Measurement of Magnetically Induced Displacement Force on Medical Devices in the Magnetic Resonance Environment
ASTM F2096-11	Standard Test Method for Detecting Gross Leaks in Packaging by Internal Pressurization
ASTM F2119-07(2013)	Standard Test Method for Evaluation of MR Image Artifacts from Passive Implants
ASTM F2182-19e2	Standard Test Method for Measurement of Radio Frequency Induced Heating On or Near Passive Implants During Magnetic Resonance Imaging

**Summary of Safety and Clinical Performance:
NuPort®/C-Port® Port Systems Revision 02**

ASTM F2213-17	Standard Test Method for Measurement of Magnetically Induced Torque on Medical Devices in the Magnetic Resonance Environment
ASTM F 3039-15	Standard Test Method for Detecting Leaks in Nonporous Packaging or Flexible Barrier Materials by Dye Penetration
MEDDEV. 2.7.1-Rev. 4 2009	Evaluation of Clinical Data: A Guide for Manufacturers and Notified Bodies
FDA Guideline	Guidance on 510(k) Submissions for Implanted Infusion Ports - 1990
GMP Brazil RDC No. 16 28.03.2013	Techn. Regulations of GMP of medical devices & in vitro diagn. devices

9. Change history

SSCP Revision number	Date	Description of the change	Revision validated by the notified body
00	13.05.2022	Initial creation of the SSCP port systems incl. accessories Author: M. Hartung, QM	☐ Validated by the Notified Body Validation language: German
01	17.03.2023	Titanium and PSU ports will be listed under two different basic UDIs from now on, therefore adjusted in the SSCP and also the article lists revised. Double lumen port and revision sets have been removed from the range and therefore from the product list. List with complications was revised and updated List of standards was updated Author: M. Meyer, QM	☐ Validated by the Notified Body Validation language: German
02	10.05.2023	List of indications and contraindications was revised and updated Basic UDI-DI of the introducer was added Table of complications was revised (general and values from studies were added) Parameters for the acceptance of safety and performance were added References to the available clinical data were added PMS and PMCF activities were supplemented Abbreviations, acronyms and explanation of medical terms were added in the patient section Library added at the end of the document Author: M. Meyer and K. Münster, QM	☒ Validated by the Notified Body Validation language: German

**Brief Safety and Clinical Performance Report:
NuPort®/C-Port® Port Systems Revision 02**

Part 2: Information for patients/public

Abbreviations and acronyms

BaSO ₄	Barium sulfate
CAPA	Corrective And Preventive Actions
CECT	Contrast Enhanced Computed Tomography
CER	Clinical Evaluation Report
CMR	Carcinogenic, Mutagenic, Reprotoxic (carcinogenic, mutagenic, fertility-endangering)
CVC	Central Venous Catheters
et al.	And others
EtO	Ethylene oxide
F	French (specification of the outer diameter of catheters)
FSCA	Field Safety Corrective Action
FSN	Field Safety Notice (safety instructions in the field)
MDR	Medical Device Regulation
ml	Milliliter
PICC	Peripherally inserted central venous catheter
PMCF	Post Marked Clinical Follow-up
PMS	Post Marked Surveillance
PSU	Polysulfone
PUR	Polyurethane
RKL	Supplier complaint
SSCP	Summary of Safety and Clinical Performance

Explanation of medical terms

Air embolism	Penetration of air into the bloodstream
Bacteremia	Presence of bacteria in the blood
Cardiac puncture	Accidental puncture of the heart
Cardiac tamponade	Accumulation of fluid in the pericardium
Catheter embolism	Closure of the catheter by a blood clot
Central venous Catheter	Thin plastic tube inserted into the venous system via a vein in the upper half of the body, the end of which lies in the superior or inferior vena cava in front of the right atrium of the heart
CE marking	In order for medical devices to be placed on the European market or put into service, they must bear a CE marking. This may only be affixed if the product meets the general safety and performance requirements.
Chemotherapy	Medicinal therapy of cancers
Chronic obstructive pulmonary disease	Disease pattern of the lungs with a permanent narrowing of the airways
Clot	Small clot of clotted blood in the bloodstream
Contraindication	Indicates when a medical treatment may not be performed
Derivative	Modified chemical compound derived from the basic substance
Embolization	Closure of a blood vessel by a blood clot detached from the site of origin.
Endocarditis	Inflammation of the inner lining of the heart

**Brief Safety and Clinical Performance Report:
NuPort®/C-Port® Port Systems Revision 02**

Endocrinology	Study of hormones, metabolism and diseases in this field
Erosion	Tissue damage on the surface
et al.	Used when a source in a scientific paper has more than two authors. The first name is given and all other names are abbreviated with et al.
Extravasation	Leakage of fluids into the tissue during the injection of drugs
Fibrin	Fibrous molecule produced by the body that closes the wound in case of injury
Fibrin sheath	Deposition of the (fibrin) molecule at the catheter outlet
Fragmentation	Fragmentation
Hematoma	Vascular injury with bleeding into the tissue (hematoma)
Indication	Indicates when medical treatment is appropriate
Infection	Penetration of pathogens into the organism
Infusion	Parenteral administration of fluids
Lymph	Fluid in which nutrients and waste are transported in the body
Migration	Migration of port or catheter in tissue or veins
Necrosis	Death of tissue
Occlusion	Occlusion of a blood vessel or the port system
Parenteral	"Bypassing the intestine"; administration of substances into the body by bypassing the intestinal tract
Penetration	Intrusion of e.g. germs into the body
Percutaneous	"Through the skin"
Perforation	Piercing or perforating vessels
Peripheral	"In the outer zones of the body [lying]"
„Pinch off"	Pinching of the catheter between the clavicle and the first rib
Pneumothorax	Penetration of air into the pleural space (between the lungs and the chest)
Port system	A subcutaneous access for infusions into the central venous system that can be used for a longer period of time
Puncture	Insertion of a cannula
Retraction	Retract
Seldinger technique	Method of puncturing blood vessels for the purpose of inserting a catheter
Sepsis	Blood poisoning
Subcutaneous	Under the skin
Superior vena cava-Syndrome	Obstruction of flow in the superior vena cava (blood vessel in the thoracic cavity, which flows into the right atrium of the heart)
Thoracic duct	Largest lymphatic vessel of the body
Thromboembolism	Disturbed blood flow due to a detached blood clot
Thrombosis	Blood clot that leads to vascular occlusion at the site of its formation

1. Labeling of the product and general information

1.1. Trade Name(s) of the Device

NuPort®/C-Port®

Table 1: Article list titanium ports

Model	Description
-------	-------------

**Brief Safety and Clinical Performance Report:
NuPort®/C-Port® Port Systems Revision 02**

CTP-XXXCP	NuPort-CT, port system set, standard, CT compatible, power injection capable, with PUR catheter
CTP-XXXIP	NuPort-CT, port system incl. introducer set, standard, CT compatible, power injection capable, with PUR catheter
CTP-XXXCS	NuPort-CT, port system set, standard, CT compatible, power injection capable, with silicone catheter
CTP-XXXIS	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
CTP-L-XXXCP	NuPort-CT, port system set, CT compatible, power injection capable, with PUR catheter, left chest placement
CTP-L-XXXIP	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with PUR catheter, left chest placement
CTP-L-XXXCS	NuPort-CT, port system set, CT compatible, power injection capable, with silicone catheter, left chest placement
CTP-L-XXXIS	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter, left chest placement
LPA-XXXCP	NuPort-LP, port system set, CT compatible, power injection capable, with PUR catheter
LPA-XXXIP	NuPort-LP, port system incl. introducer set, CT compatible, power injection capable, with PUR catheter
LPA-XXXCS	NuPort-LP, port system set, CT compatible, power injection capable, with silicone catheter
LPA-XXXIS	NuPort-LP, port system incl. introducer set, CT-compatible, power injection capable, with silicone catheter
SPB-XXXCP	C-Port, port system set, standard, with PUR catheter
SPB-XXXIP	C-Port, port system incl. introducer set, standard, with PUR catheter
SPB-XXXCS	C-Port, port system set, standard, with silicone catheter
SPB-XXXIS	C-Port, port system incl. introducer set, standard, with silicone catheter
PH-XXIS-YY	NuPort-CT, port system incl. introducer set, standard, with silicone catheter
LPA-XXXIP-YY	NuPort-LPA, port system incl. introducer set, with PUR catheter
CTP-XXXCP-YYY	NuPort-CT, port system set, CT compatible, with PUR catheter
CTP-XXXIP-YYY	NuPort-CT, port system incl. introducer set, CT compatible, with PUR catheter
CTP-L-XXXCP-YYY	NuPort-CT, port system set, CT compatible, with PUR catheter, left chest placement
CTP-L-XXXIP-YYY	NuPort-CT, port system incl. introducer set, CT compatible, with PUR catheter, left chest placement
8086551	Fresenius C-Port HF 8F, 8F silicone catheter, with introducer set
8086561	Fresenius C-Port HF 9.6F, 9.6F silicone catheter, with introducer set

Table 2: Article list PSU ports

Model	Description
CTKP-XXXCP	C-Port-CT, port system set, standard, CT compatible, power injection capable, with PUR catheter
CTKP-XXXIP	C-Port-CT, port system incl. introducer set, standard, CT compatible, power injection capable, with PUR catheter

Brief Safety and Clinical Performance Report: NuPort®/C-Port® Port Systems Revision 02

CTKP-XXXCS	C-Port-CT, port system set, standard, CT compatible, power injection capable, with silicone catheter
CTKP-XXXIS	C-Port-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
KP-XXXCP	C-Port-II, port system set, standard, with PUR catheter
KP-XXXIP	C-Port-II, port system incl. introducer set, standard, with PUR catheter
KP-XXXCS	C-Port-II, port system set, standard, with silicone catheter
KP-XXXIS	C-Port-II, port system incl. introducer set, standard, with silicone catheter
CTKP-XXXIS-YY	C-Port-CT, port system incl. introducer set, standard, with silicone catheter
CTKP-XXXIP-YY	C-Port-CT, port system incl. introducer set, standard, with PUR catheter
CTKP-XXXCP-YYY	C-Port-CT, port system set, CT compatible, with PUR catheter

XXX - catheter size

YY/Y - abbreviation for customers of private label port systems (distribution by customers, placing on the market by PHS Medical GmbH)

The ports are available with catheters in different sizes. The following catheters are available:

- 5F PUR
- 6.6F PUR
- 7.5F PUR
- 9F PUR
- 8F Silicone
- 9.6F Silicone

1.2. Name and address of the manufacturer

PHS Medical GmbH
Ederweg 3
34277 Fuldabrück
Germany

Tel.: +49 (0) 561 - 99 85 97 0

E-mail: info@phs-medical.de

1.3. Basic UDI-DI

Titanium ports: 42503746PortsystemTitanNF

PSU ports: 42503746PortsystemPSUC6

1.4. Year in which the first certificate (CE) for the product was issued

2003

2. Purpose of the product

**Brief Safety and Clinical Performance Report:
NuPort®/C-Port® Port Systems Revision 02**

2.1. Intended use

The **NuPort®/C-Port®** port family are medical devices that are implanted under the skin, mostly in the chest area, to provide access to the vascular system. These systems are designed for patients who need to receive repeated or regular medication through the vein.

These products provide medical personnel with an easy-to-find site for inserting a cannula to administer fluids or drugs into the vasculature or to collect blood samples. The port is accessed by puncture (piercing) with a (non-coring) port cannula.

Ports are intended for single use and are explanted and disposed of at the end of therapy.

2.2. Indications and target group

The **NuPort®/C-Port®** family of ports are intended for use with patients who require the following:

- repeated access to the vascular system for injections,
- infusion of drugs,
- administration of blood or blood products and/or collection of blood as part of therapy.

Areas of application include:

- Parenteral nutrition
- Chemotherapy
- Pain therapy
- endocrinological therapy
- Substitution of coagulation factors
- Contrast medium infusion in radiation diagnostics
- Therapy of severe asthma or cardiac arrhythmias
- Blood collection and transfusion

NuPort®/C-Port® are intended for all patients who require long-term vascular access. The physician responsible for you will select the port model according to your age and body stature.

2.3. Contraindication and restrictions

The ports of the **NuPort®/C-Port®** family are not suitable for use under the following conditions:

- Patients with known or suspected allergic reactions to the materials contained in the implant
- Patients with known infection, bacteremia, or sepsis
- Patients with body tissue that cannot adequately support the port or catheter
- Patients with chronic obstructive pulmonary disease (COPD)
- When a potential access site has been irradiated
- If the patient's tissue factors prevent proper stabilization of the product and/or access
- Previous venous thrombosis or surgical procedures at or near the site of the planned use

3. Product description

3.1. Description of the medical device

**Brief Safety and Clinical Performance Report:
NuPort®/C-Port® Port Systems Revision 02**

The product is a set of various components which are sterilized using EtO (ethylene oxide). It is a single-use product, which must not be reused.

The **implantable** parts include a port, which is made of titanium or polysulfone (PSU) depending on the variant, and in which a silicone septum is mounted, as well as a catheter made of silicone or polyurethane (PUR) and a catheter fixation, which are made of polysulfone or polycarbonate. The septum is manufactured to withstand a very high number of punctures with a non-coring cannula. To connect the catheter to the port, a connector is screwed into the port and glued with silicone adhesive. The catheter is pushed onto this during implantation and fixed with the catheter fixation for a secure connection.

For implantation, the following components for **single use** are also included in the set:

- Huber cannula, which enables a puncture of the septum without coring
(Material: cannula - stainless steel, luer - polycarbonate)
- Blunt cannula for flushing the catheter
(Material: cannula - stainless steel, luer - polycarbonate)
- 10ml syringe with luer slip
(Material: polypropylene, polyisoprene)
- A vein pick, which is designed to facilitate the insertion of the catheter into the vein
(Material: terlux)

Furthermore, the portsets are also available as a variant with introducer. This is used for insertion of the catheters via the minimally invasive percutaneous technique (cf. "Seldinger technique") or in the so-called open surgical technique.

The following additional components are included in the set with Introducer:

- Puncture cannula (Seldinger cannula)
(Material: cannula - stainless steel, luer - polycarbonate)
- Guide wire with dispenser and insertion aid
(Material: stainless steel, polypropylene, polyethylene)
- Split sheath (tearaway) with dilator
(Material: polyflon, polypropylene, polyethylene)
- Tunneler
(Material: Stainless steel)

**Brief Safety and Clinical Performance Report:
NuPort®/C-Port® Port Systems Revision 02**


Fig 1: Portset with Titanium port and Introducer

3.2. Information on the medicinal products contained in the product, if any

The ports of the **NuPort®/C-Port®** family do **not** contain **any**

- Drug substances (including derivatives from human blood or plasma)
- Tissues or cells of human or animal origin or their derivatives
- Substances or combinations of substances that are absorbed by the human body or distributed locally in the body
- Materials containing CMR substances (carcinogens, mutagens or substances toxic to reproduction) or endocrine disruptors.

While silicone is a known allergen, there is no evidence of allergies to medical grade silicone in the literature or from practice.

3.3. Description of the way in which the product achieves its intended effect

The port is a device approximately 3 cm in size that is inserted completely under the skin in a short procedure. The port system was developed for patients who need to receive medication repeatedly or regularly via a vein. The base consists of a cylinder with a cavity (reservoir), which is closed at the top with an elastic membrane (septum). A small, flexible tube (catheter) is connected to the port and inserted into a large central vein that carries blood to the heart. The system is implanted completely under the skin and provides secure access to a large vein that can be used at any time. A special, non-coring port cannula, which is inserted into the septum of the port, is used to "access" the bloodstream. Drugs and fluids can be administered and blood samples taken in this way.

3.4. Description of the accessories, if any

**Brief Safety and Clinical Performance Report:
NuPort[®]/C-Port[®] Port Systems Revision 02**

There are no accessories relevant to the patients.

4. Risks and undesirable effects

Consult your healthcare professional if you think you may experience side effects related to the product or its use, or if you are concerned about any risks. This document is not a substitute for advice from a healthcare professional.

4.1. How were potential risks controlled or managed?

Potential risks were identified and controlled in accordance with the standard DIN EN ISO 14971 Medical devices - Application of risk management to medical devices.

4.2. Residual risks and undesirable effects

Use of the **NuPort[®]/C-Port[®]** family of ports involves risks normally associated with the insertion or use of an implant or indwelling catheter, including but not limited to the following:

- Infection
- Erosion of the skin
- Hematoma
- Necrosis
- Displacement of the product
- Clot formation
- Catheter fragmentation
- Embolization
- Catheter closure
- Cardiac arrhythmias
- Cardiac puncture
- Pericardial tamponade
- Extravasation of drugs
- Fibrin sheath
- Endocarditis
- Implant rejections
- Port/Catheter Migration
- Perforation or rupture of vessels
- Injury of the thoracic duct
- Thromboembolism
- Thrombosis
- Superior vena cava syndrome
- Pneumothorax
- spontaneous malpositioning or retraction of the catheter tip
- Vascular erosion
- catheter- or port-related sepsis
- Air embolism

**Brief Safety and Clinical Performance Report:
NuPort®/C-Port® Port Systems Revision 02**

- Bleeding
- Catheter embolism

The following table shows the most frequent complications reported to PHS Medical GmbH in the last 20 years and their frequency in scientific literature and from studies performed on port systems of PHS Medical GmbH:

Table 3: Frequency of complications of port systems

Complication	Average Complication rates from the literature*	Average complication rates from studies performed on products of PHS Medical GmbH	Number of complications reported to PHS Medical GmbH (2002-2022)	Percentage share in relation to the total volume of products sold to date at PHS Medical GmbH (2002- 2022)
Complications due to disregard of the instructions for use	Not named	0,5%	39	0,011%
"Pinch off", catheter rupture and dislocation.	1%	0,9%	29	0,013%
Damaged septum	Not named	0%	27	0,012%
Occlusion	0,9%	0,9%	2	0,001%
Infection	4,5%	4,2%	1	0,0003%
Thromboses	2,8%	0,9%	0	0%

4.3. Warnings and precautions

Patient Warnings and Precautions:

- Inform your doctor if you have any allergies to medications or other substances so that the doctor can select the port or medications to use accordingly.
- Inform your doctor if you have any pre-existing conditions or have received radiation therapy.
- Avoid physical exertion for the first few days after implantation
- Follow the instructions of your doctor or nurse for the care of the wound suture
- Avoid activities that may cause bruising in the area of the port. This may result in damage to the catheter.
- Please inform your doctor immediately if:
 - You notice redness or inflammation at the site of the port after the incision has healed, or
 - You get a fever.
- Consult your doctor if you want to do any special activities, such as diving, after healing. You can resume everyday activities such as swimming or jogging.
- A bandage is only required until the incision healed, after that it is no longer necessary.
- The port may only be punctured by medical personnel.
- Inform the medical staff if you experience pain during the infusion.

**Brief Safety and Clinical Performance Report:
NuPort®/C-Port® Port Systems Revision 02**

- During breaks in therapy, PHS Medical GmbH recommends flushing the port system with saline solution every three months.

Various warnings and precautions must be observed during implantation and use. A detailed description for the users (medical professionals) can be found in Part 1 of this document.

Note

The user is obliged to report all serious incidents related to the product to the manufacturer and the competent authorities in the country where the user is located.

4.4. Other relevant safety aspects, including a summary of all field safety corrective actions (FSCA including FSN), if applicable.

In 2017, a recall occurred due to the lack of X-ray visibility of the F9.6 silicone catheters. This incident was registered at the Federal Institute for Medical Devices and Drugs under case number 01229/17, and at PHS Medical GmbH under the internal complaint number RKL036 and in CAPA002. The processing of the incident was completed on 07 April 2017. All products containing affected catheters have been recalled.

The cause of the event was insufficient mixing of the materials. Although the silicone base material and BaSO₄ were correctly weighed and measured, the mixture - according to the operator - was not thoroughly processed with the intended 25 cycles in the grinder. Therefore, the BaSO₄ distribution was not homogeneous.

A CAPA was opened at the catheter supplier and the relevant work instructions were adapted to ensure that the materials were sufficiently mixed together.

In the meantime, the supplier for the silicone catheters has been changed.

5. Summary of clinical evaluation and post-market clinical follow-up

5.1. Clinical background of the product

Central venous port catheter systems are becoming increasingly important today, especially in oncology, as they provide permanent vascular access, significantly facilitate the administration of drugs, chemotherapeutics, high-caloric nutrient solutions and blood products (etc.) as well as blood sampling, and minimize the risk of infection due to their location under the skin. High-pressure-tested port systems can also be used as part of a contrast-enhanced CT scan (power injected contrast injection).

The first long-term venous catheters were available in 1968 and were made of polyethylene at that time. Then, in the early 1970s, Broviac et al. described for the first time the use of softer and more flexible silicone for central venous catheters for parenteral nutrition therapy. This was the first significant advance in long-term catheters. In the late 1970s, Hickman et al reported modifications, including a larger catheter lumen. However, the catheter segments protruding from the skin remained a disadvantage because they provided an access route for germs and an associated high risk of infection from skin penetration. The exit site of these catheter variants, which still requires care, and the limited freedom of movement have also remained to this day. Then, in 1982, Niederhuber et al described the clinical introduction of subcutaneously (under the skin) implanted port systems for chemotherapy. Based on these encouraging results, port catheter systems became more widespread in the following years, and reports of experience with the use of such systems, including in larger patient populations, increased.

Brief Safety and Clinical Performance Report: NuPort®/C-Port® Port Systems Revision 02

The implantable ports available today can be considered technically mature. In this context, the ports available on the market are increasingly similar in shape and components.

Most port manufacturers, such as PHS Medical GmbH, offer different ports whose size and height are adapted to the respective purpose (ports for adults, mini-ports for children and underweight patients). This allows the port to be found during puncture without being unnecessarily bulky. The large septum (at least as large as a fingertip) and the palpable contour of the port (rim bulge around the membrane) make it easier to perform the port puncture. The localization of the port can be detected and thus a safe injection is ensured during use.

The shape of the port should be such that the port can be easily palpated. The shape should not create unnecessary subcutaneous cavities. Sharp edges should be avoided because of the risk of skin erosion. The interior of the port should be smooth and have no corners where blood components can settle. This plays a role especially with regard to the prevention of infections inside the port. The more sharp edges there are inside the port, the greater the risk of blood components being deposited, which in turn provide a good breeding ground for germs of all kinds. The port systems from PHS Medical GmbH are designed to meet all these points.

One parameter for demonstrating the performance of PHS Medical GmbH products is their ability to provide such long-term vascular access ("Long-term" > 30 days). Here, a permanently functioning seal between the port chamber and tissue and the transfer of fluids between the port chamber and the outlet site of the catheter should prevent irritation/damage to surrounding tissue and veins by the administered drugs. The lifetime of port systems is limited by the number of punctures with a cannula (1000/cm²).

The reported complication rates from the literature are used as a basis for evaluating the clinical safety of the products from PHS Medical GmbH. Depending on the source, these range between 2.6 and 46 %. Taking into account the average number of complaints about the port systems and introducer sets from PHS Medical GmbH of less than 1%, an overall internal complication rate of max. 2% for the port systems incl. accessories and introducer sets was therefore determined or defined as acceptable.

Both the clinical studies on the port systems and introducer sets from PHS Medical GmbH and the overall assessment of the port systems and introducer sets from PHS Medical GmbH in the clinical evaluation show that the **NuPort®/C-Port®** port systems from PHS Medical GmbH under consideration provide safe and long-term access to a patient's vascular system for a variety of indications. The established parameters for demonstrating safety and performance are met by the PHS Medical GmbH products.

The development and manufacture of the products are carried out in accordance with the current legal requirements and thus meet the safety and performance requirements.

According to current knowledge and the results of tests, clinical data and scientific literature, the benefits outweigh the risks of using the port systems and introducer sets. The analysis and evaluation of potential risks has shown that there are no increased residual risks for patients, users or third parties within the scope of the intended use of the port systems and introducer sets.

5.2. The clinical evidence for CE marking

The clinical evidence for the CE marking of the **NuPort®/C-Port®** port systems is based on laboratory tests, scientific literature, market feedback and own clinical data on the products from clinical studies. Studies

Brief Safety and Clinical Performance Report: NuPort®/C-Port® Port Systems Revision 02

conducted on PHS Medical GmbH products at the Universities of Münster and Essen, among others, demonstrated that the **NuPort®/C-Port®** port systems function according to their intended use and are safe for use on patients. In both studies, it was demonstrated on a large number of patients that the implantation and use of port systems is a safe procedure to obtain permanent venous access.

In the study of the University of Essen by Schenck et al. (2012), 347 patients were implanted with **NuPort®** (titanium port with silicone catheter) from PHS Medical GmbH between December 1999 and June 2011. In most cases (91.7%), implantation was due to planned chemotherapy. The other applications were artificial nutrition (4.6%) and lack of suitable peripheral vascular access (1.4%) or pain therapy (2.3%). The implantations of the port systems could be performed without complications in all cases. Early explantation of the port system within the first 30 days was necessary in six cases, of which three were due to catheter migration and one each to sepsis, infection, and thrombosis. The average duration of stay was 239 days (2 - 2604 days), demonstrating that the port systems from PHS Medical GmbH are suitable to provide permanent (> 30 days) access to the vascular system.

This is also confirmed by the study of the University of Münster by Becker et al. (2021), who implanted a port system from PHS Medical GmbH in 500 patients and followed the patients over a period of 12 months. Introducer sets from PHS Medical GmbH were used for implantation. Implantations were mostly complication-free (4 complications 0.8%). In one case, the goal of establishing long-term (> 30 days) access to the vascular system was not achieved because the port system had to be explanted prematurely due to an infection.

Overall, the results presented in these studies were comparable to those in the literature in terms of success and complications, demonstrating that the **NuPort®/C-Port®** port systems perform according to their intended use and are safe to use.

The data from these studies are incorporated into the clinical evaluation of the safety and performance of the products from PHS Medical GmbH. This clinical evaluation shows that the **NuPort®/C-Port®** port systems under consideration provide safe and long-term access to a patient's vascular system for a variety of indications.

According to the current state of knowledge, based on the state of the art as well as the product-specific results from tests, clinical data and scientific literature, the benefits outweigh the risks of using the port systems. The analysis and evaluation of potential risks has shown that there are no increased residual risks for patients, users or third parties in the context of the intended use of the port systems and introducer sets, which can be confirmed by the product-related clinical data.

The benefit of the **NuPort®/C-Port®** port systems and introducer sets for the patients has been clearly demonstrated.

5.3. Safety

Port systems have been used in humans for many years and the complications and side effects associated with their use are well documented.

The clinical evaluation of PHS Medical GmbH shows that the **NuPort®/C-Port®** port systems are safe and effective and do not endanger the clinical condition and safety of patients or the safety and health of users or third parties, if applicable.

The safety and performance of PHS Medical GmbH products are continuously and actively monitored. For this purpose, data is collected on a regular basis. These so-called PMS (market surveillance) and PMCF (post market clinical follow-up) activities, which are carried out at PHS Medical GmbH, include:

- Customer surveys on satisfaction with port systems as well as product- and safety-related topics
- Review and evaluation of complaints and CAPAs

**Brief Safety and Clinical Performance Report:
NuPort[®]/C-Port[®] Port Systems Revision 02**

- Review and evaluation of observations, reports, events and occurrences made to PHS Medical GmbH.
- Audit of risk management
- Inspection and evaluation of data from the production process
- Review of the instructions for use
- Review of relevant new and revised standards
- Review of the relevant scientific literature
- Checking relevant databases and registers
- Recording of all data in the PHS Medical GmbH trend analysis, regular review of this data and the statistics contained therein
- Review of the PMCF plan and report.

The results from this market monitoring are recorded in an annual safety report. Due to the long market presence of the **NuPort[®]/C-Port[®]** port systems as well as the clinical data already available, no clinical studies are currently being conducted as part of the PMCF activities.

The evaluation performed in the current safety report and PMCF report has confirmed the safety and performance of the PHS port catheter systems during their expected lifetime. No previously unknown adverse reactions, side effects, contraindications, or emerging risks were identified during the surveillance period. It was confirmed that the benefit-risk ratio is acceptable.

6. Therapeutic alternatives

Only central venous catheters (CVC) and peripherally inserted central venous catheters (PICC) represent an alternative to port systems. However, port systems have the advantage over central venous catheters that they are less susceptible to infection, require less care, can be used for longer periods of time, have higher cosmetic acceptability, and are perceived as less burdensome by the patient. To prevent infection with a CVC, the insertion site must be cleaned daily. In addition, complications such as catheter rupture or thrombosis can also occur with the CVC.

If you are considering alternative treatments, you should consult your doctor, who can take your individual situation into account.

7. User qualification and training

The **NuPort[®]/C-Port[®]** port systems may only be used by a physician or medical staff. Therefore, training for patients is not necessary.

**Brief Safety and Clinical Performance Report:
NuPort[®]/C-Port[®] Port Systems Revision 02****Appendix I: Library:**

Becker, F., Wurche, L. A., Darscht, M., Pascher, A., & Struecker, B. (2021). Totally implantable venous access port insertion via open Seldinger approach of the internal jugular vein-a retrospective risk stratification of 500 consecutive patients. *Langenbeck's Archives of Surgery*, 406(3), 903-910.

Busch, Jasmin D., et. al, Complication Rates Observed in Silicone and Polyurethane Catheters of TIVADs Implanted in the Upper Arm, *J Vasc Interv Radiol* 2017, 1-7.

Fenner, Katharina, Duration of use and complications of intravenous port systems - a retrospective study of 465 ports of the Clinic and Polyclinic for Anaesthesiology and Intensive Therapy at the University of Leipzig from 1998 - 2001, Medical Faculty of the University of Leipzig 2005: Dissertation.

Jordan, Karin, et. A., Venous Access Ports: Frequency and Management of Complications in Oncology Patients, *Oncology* 2008;31:404-410.

Kalmbach, Sarah Johanna, Risk factors for premature port catheter explantation - a complication analysis of more than 1 000 port catheters, Eberhard Karls University of Tübingen 2021: Dissertation.

Kunz, Julia Rebecca, Evaluation of success and complication rates of radiologically implanted power injection port systems depending on vascular access and experience of the interventional radiologist, Faculty of Medicine Charité - Universitätsmedizin Berlin 2018: Dissertation.

Kunz-Virk, Julia, Krüger, Karsten, Power-injection totally implantable venous access devices - analysis of success and complication rates of ultrasound-guided implantation, *Vasa* (2019), 48, 524-530.

Lippert, Hans, Mantke, René (2018), Risks and complications in general and visceral surgery, Stuttgart; Thieme Verlag 2018, p. 204.

Schenck M, Michels-Oswald W, Tschirdewahn S, Rübber H, Dorp F vom, Rose A, et al. How should urologists implant venous subcutaneous port systems? *Urol. Feb.* 1, 2012;51(2):226-37.

Völklein, Jana, Minimally Invasive Central Venous Port Implantations - Comparison of Complications of Surgical Technique, Ludwig Maximilian University, Munich 2020: Dissertation.

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Modelo	Descripción
CTP-XXXCP	NuPort-CT, port system set, standard, CT compatible, power injection capable, with PUR catheter
CTP-XXXIP	NuPort-CT, port system incl. introducer set, standard, CT compatible, power injection capable, with PUR catheter
CTP-XXXCS	NuPort-CT, port system set, standard, CT compatible, power injection capable, with silicone catheter
CTP-XXXIS	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
CTP-L-XXXCP	NuPort-CT, port system set, CT compatible, power injection capable, with PUR catheter, left chest placement
CTP-L-XXXIP	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with PUR catheter, left chest placement
CTP-L-XXXCS	NuPort-CT, port system set, CT compatible, power injection capable, with silicone catheter, left chest placement
CTP-L-XXXIS	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter, left chest placement
LPA-XXXCP	NuPort-LP, low-profile port system set, CT compatible, power injection capable, with PUR catheter
LPA-XXXIP	NuPort-LP, low-profile port system incl. introducer set, CT compatible, power injection capable, with PUR catheter
LPA-XXXCS	NuPort-LP, low-profile port system set, CT compatible, power injection capable, with silicone catheter
LPA-XXXIS	NuPort-LP, low-profile port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
SPB-XXXCP	C-Port, port system set, standard, with PUR catheter
SPB-XXXIP	C-Port, port system incl. introducer set, standard, with PUR catheter
SPB-XXXCS	C-Port, port system set, standard, with silicone catheter
SPB-XXXIS	C-Port, port system incl. introducer set, standard, with silicone catheter
PH-XXIS-YY	NuPort-CT, port system incl. introducer set, standard, with silicone catheter
LPA-XXXIP-YY	NuPort-LPA, low-profile port system incl. introducer set, with PUR catheter
CTP-XXXCP-YYY	NuPort-CT, port system set, CT compatible, with PUR catheter
CTP-XXXIP-YYY	NuPort-CT, port system incl. introducer set, CT compatible, with PUR catheter

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CTP-L-XXXCP-YYY	NuPort-CT, port system set, CT compatible, with PUR catheter, left chest placement
CTP-L-XXXIP-YYY	NuPort-CT, port system incl. introducer set, CT compatible, with PUR catheter, left chest placement
8086551	Fresenius C-Port HF 8F, 8F silicone catheter, with introducer set
8086561	Fresenius C-Port HF 9.6F, 9.6F silicone catheter, with introducer set

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Modelo	Descripción
CTKP-XXXCP	C-Port-CT, port system set, standard, CT compatible, power injection capable, with PUR catheter
CTKP-XXXIP	C-Port-CT, port system incl. introducer set, standard, CT compatible, power injection capable, with PUR catheter
CTKP-XXXCS	C-Port-CT, port system set, standard, CT compatible, power injection capable, with silicone catheter
CTKP-XXXIS	C-Port-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
KP-XXXCP	C-Port-II, port system set, standard, with PUR catheter
KP-XXXIP	C-Port-II, port system incl. introducer set, standard, with PUR catheter
KP-XXXCS	C-Port-II, port system set, standard, with silicone catheter
KP-XXXIS	C-Port-II, port system incl. introducer set, standard, with silicone catheter
CTKP-XXXIS-YY	C-Port-CT, port system incl. introducer set, standard, with silicone catheter
CTKP-XXXIP-YY	C-Port-CT, port system incl. introducer set, standard, with PUR catheter
CTKP-XXXCP-YYY	C-Port-CT, port system set, CT compatible, with PUR catheter

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Modelo	Descripción
CTP-XXXCP	NuPort-CT, port system set, standard, CT compatible, power injection capable, with PUR catheter
CTP-XXXIP	NuPort-CT, port system incl. introducer set, standard, CT compatible, power injection capable, with PUR catheter
CTP-XXXCS	NuPort-CT, port system set, standard, CT compatible, power injection capable, with silicone catheter
CTP-XXXIS	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
CTP-L-XXXCP	NuPort-CT, port system set, CT compatible, power injection capable, with PUR catheter, left chest placement
CTP-L-XXXIP	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with PUR catheter, left chest placement
CTP-L-XXXCS	NuPort-CT, port system set, CT compatible, power injection capable, with silicone catheter, left chest placement
CTP-L-XXXIS	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter, left chest placement
LPA-XXXCP	NuPort-LP, port system set, CT compatible, power injection capable, with PUR catheter
LPA-XXXIP	NuPort-LP, port system incl. introducer set, CT compatible, power injection capable, with PUR catheter
LPA-XXXCS	NuPort-LP, port system set, CT compatible, power injection capable, with silicone catheter
LPA-XXXIS	NuPort-LP, port system incl. introducer set, CT-compatible, power injection capable, with silicone catheter
SPB-XXXCP	C-Port, port system set, standard, with PUR catheter
SPB-XXXIP	C-Port, port system incl. introducer set, standard, with PUR catheter
SPB-XXXCS	C-Port, port system set, standard, with silicone catheter
SPB-XXXIS	C-Port, port system incl. introducer set, standard, with silicone catheter
PH-XXIS-YY	NuPort-CT, port system incl. introducer set, standard, with silicone catheter
LPA-XXXIP-YY	NuPort-LPA, port system incl. introducer set, with PUR catheter
CTP-XXXCP-YYY	NuPort-CT, port system set, CT compatible, with PUR catheter
CTP-XXXIP-YYY	NuPort-CT, port system incl. introducer set, CT compatible, with PUR catheter
CTP-L-XXXCP-YYY	NuPort-CT, port system set, CT compatible, with PUR catheter, left chest placement
CTP-L-XXXIP-YYY	NuPort-CT, port system incl. introducer set, CT compatible, with PUR catheter, left chest placement
8086551	Fresenius C-Port HF 8F, 8F silicone catheter, with introducer set
8086561	Fresenius C-Port HF 9.6F, 9.6F silicone catheter, with introducer set

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CTKP-XXXCP	C-Port-CT, port system set, standard, CT compatible, power injection capable, with PUR catheter
CTKP-XXXIP	C-Port-CT, port system incl. introducer set, standard, CT compatible, power injection capable, with PUR catheter
CTKP-XXXCS	C-Port-CT, port system set, standard, CT compatible, power injection capable, with silicone catheter
CTKP-XXXIS	C-Port-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
KP-XXXCP	C-Port-II, port system set, standard, with PUR catheter
KP-XXXIP	C-Port-II, port system incl. introducer set, standard, with PUR catheter
KP-XXXCS	C-Port-II, port system set, standard, with silicone catheter
KP-XXXIS	C-Port-II, port system incl. introducer set, standard, with silicone catheter
CTKP-XXXIS-YY	C-Port-CT, port system incl. introducer set, standard, with silicone catheter
CTKP-XXXIP-YY	C-Port-CT, port system incl. introducer set, standard, with PUR catheter
CTKP-XXXCP-YYY	C-Port-CT, port system set, CT compatible, with PUR catheter

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Μοντέλο	Περιγραφή
CTP-XXXCP	NuPort-CT, port system set, standard, CT compatible, power injection capable, with PUR catheter
CTP-XXXIP	NuPort-CT, port system incl. introducer set, standard, CT compatible, power injection capable, with PUR catheter
CTP-XXXCS	NuPort-CT, port system set, standard, CT compatible, power injection capable, with silicone catheter
CTP-XXXIS	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
CTP-L-XXXCP	NuPort-CT, port system set, CT compatible, power injection capable, with PUR catheter, left chest placement
CTP-L-XXXIP	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with PUR catheter, left chest placement
CTP-L-XXXCS	NuPort-CT, port system set, CT compatible, power injection capable, with silicone catheter, left chest placement
CTP-L-XXXIS	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter, left chest placement
LPA-XXXCP	NuPort-LP, low-profile port system set, CT compatible, power injection capable, with PUR catheter
LPA-XXXIP	NuPort-LP, low-profile port system incl. introducer set, CT compatible, power injection capable, with PUR catheter
LPA-XXXCS	NuPort-LP, low-profile port system set, CT compatible, power injection capable, with silicone catheter
LPA-XXXIS	NuPort-LP, low-profile port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
SPB-XXXCP	C-Port, port system set, standard, with PUR catheter
SPB-XXXIP	C-Port, port system incl. introducer set, standard, with PUR catheter
SPB-XXXCS	C-Port, port system set, standard, with silicone catheter
SPB-XXXIS	C-Port, port system incl. introducer set, standard, with silicone catheter
PH-XXIS-YY	NuPort-CT, port system incl. introducer set, standard, with silicone catheter
LPA-XXXIP-YY	NuPort-LPA, low-profile port system incl. introducer set, with PUR catheter
CTP-XXXCP-YYY	NuPort-CT, port system set, CT compatible, with PUR catheter
CTP-XXXIP-YYY	NuPort-CT, port system incl. introducer set, CT compatible, with PUR catheter

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Μοντέλο	Περιγραφή
CTKP-XXXCP	C-Port-CT, port system set, standard, CT compatible, power injection capable, with PUR catheter
CTKP-XXXIP	C-Port-CT, port system incl. introducer set, standard, CT compatible, power injection capable, with PUR catheter
CTKP-XXXCS	C-Port-CT, port system set, standard, CT compatible, power injection capable, with silicone catheter
CTKP-XXXIS	C-Port-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
KP-XXXCP	C-Port-II, port system set, standard, with PUR catheter
KP-XXXIP	C-Port-II, port system incl. introducer set, standard, with PUR catheter
KP-XXXCS	C-Port-II, port system set, standard, with silicone catheter
KP-XXXIS	C-Port-II, port system incl. introducer set, standard, with silicone catheter
CTKP-XXXIS-YY	C-Port-CT, port system incl. introducer set, standard, with silicone catheter
CTKP-XXXIP-YY	C-Port-CT, port system incl. introducer set, standard, with PUR catheter
CTKP-XXXCP-YYY	C-Port-CT, port system set, CT compatible, with PUR catheter

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Μοντέλο	Περιγραφή
CTP-XXXCP	NuPort-CT, port system set, standard, CT compatible, power injection capable, with PUR catheter
CTP-XXXIP	NuPort-CT, port system incl. introducer set, standard, CT compatible, power injection capable, with PUR catheter
CTP-XXXCS	NuPort-CT, port system set, standard, CT compatible, power injection capable, with silicone catheter
CTP-XXXIS	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
CTP-L-XXXCP	NuPort-CT, port system set, CT compatible, power injection capable, with PUR catheter, left chest placement
CTP-L-XXXIP	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with PUR catheter, left chest placement
CTP-L-XXXCS	NuPort-CT, port system set, CT compatible, power injection capable, with silicone catheter, left chest placement
CTP-L-XXXIS	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter, left chest placement
LPA-XXXCP	NuPort-LP, port system set, CT compatible, power injection capable, with PUR catheter
LPA-XXXIP	NuPort-LP, port system incl. introducer set, CT compatible, power injection capable, with PUR catheter
LPA-XXXCS	NuPort-LP, port system set, CT compatible, power injection capable, with silicone catheter
LPA-XXXIS	NuPort-LP, port system incl. introducer set, CT-compatible, power injection capable, with silicone catheter
SPB-XXXCP	C-Port, port system set, standard, with PUR catheter
SPB-XXXIP	C-Port, port system incl. introducer set, standard, with PUR catheter
SPB-XXXCS	C-Port, port system set, standard, with silicone catheter
SPB-XXXIS	C-Port, port system incl. introducer set, standard, with silicone catheter
PH-XXIS-YY	NuPort-CT, port system incl. introducer set, standard, with silicone catheter
LPA-XXXIP-YY	NuPort-LPA, port system incl. introducer set, with PUR catheter

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CTP-XXXCP-YYY	NuPort-CT, port system set, CT compatible, with PUR catheter
CTP-XXXIP-YYY	NuPort-CT, port system incl. introducer set, CT compatible, with PUR catheter
CTP-L-XXXCP-YYY	NuPort-CT, port system set, CT compatible, with PUR catheter, left chest placement
CTP-L-XXXIP-YYY	NuPort-CT, port system incl. introducer set, CT compatible, with PUR catheter, left chest placement
8086551	Fresenius C-Port HF 8F, 8F silicone catheter, with introducer set
8086561	Fresenius C-Port HF 9.6F, 9.6F silicone catheter, with introducer set

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Modelle	Descrizione
CTP-XXXCP	NuPort-CT, port system set, standard, CT compatible, power injection capable, with PUR catheter
CTP-XXXIP	NuPort-CT, port system incl. introducer set, standard, CT compatible, power injection capable, with PUR catheter
CTP-XXXCS	NuPort-CT, port system set, standard, CT compatible, power injection capable, with silicone catheter
CTP-XXXIS	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
CTP-L-XXXCP	NuPort-CT, port system set, CT compatible, power injection capable, with PUR catheter, left chest placement
CTP-L-XXXIP	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with PUR catheter, left chest placement
CTP-L-XXXCS	NuPort-CT, port system set, CT compatible, power injection capable, with silicone catheter, left chest placement
CTP-L-XXXIS	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter, left chest placement
LPA-XXXCP	NuPort-LP, low-profile port system set, CT compatible, power injection capable, with PUR catheter
LPA-XXXIP	NuPort-LP, low-profile port system incl. introducer set, CT compatible, power injection capable, with PUR catheter
LPA-XXXCS	NuPort-LP, low-profile port system set, CT compatible, power injection capable, with silicone catheter
LPA-XXXIS	NuPort-LP, low-profile port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
SPB-XXXCP	C-Port, port system set, standard, with PUR catheter
SPB-XXXIP	C-Port, port system incl. introducer set, standard, with PUR catheter
SPB-XXXCS	C-Port, port system set, standard, with silicone catheter
SPB-XXXIS	C-Port, port system incl. introducer set, standard, with silicone catheter
PH-XXIS-YY	NuPort-CT, port system incl. introducer set, standard, with silicone catheter
LPA-XXXIP-YY	NuPort-LPA, low-profile port system incl. introducer set, with PUR catheter
CTP-XXXCP-YYY	NuPort-CT, port system set, CT compatible, with PUR catheter
CTP-XXXIP-YYY	NuPort-CT, port system incl. introducer set, CT compatible, with PUR catheter

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Modelle	Descrizione
CTKP-XXXCP	C-Port-CT, port system set, standard, CT compatible, power injection capable, with PUR catheter
CTKP-XXXIP	C-Port-CT, port system incl. introducer set, standard, CT compatible, power injection capable, with PUR catheter
CTKP-XXXCS	C-Port-CT, port system set, standard, CT compatible, power injection capable, with silicone catheter
CTKP-XXXIS	C-Port-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
KP-XXXCP	C-Port-II, port system set, standard, with PUR catheter
KP-XXXIP	C-Port-II, port system incl. introducer set, standard, with PUR catheter
KP-XXXCS	C-Port-II, port system set, standard, with silicone catheter
KP-XXXIS	C-Port-II, port system incl. introducer set, standard, with silicone catheter
CTKP-XXXIS-YY	C-Port-CT, port system incl. introducer set, standard, with silicone catheter
CTKP-XXXIP-YY	C-Port-CT, port system incl. introducer set, standard, with PUR catheter
CTKP-XXXCP-YYY	C-Port-CT, port system set, CT compatible, with PUR catheter

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Modelle	Descrizione
CTP-XXXCP	NuPort-CT, port system set, standard, CT compatible, power injection capable, with PUR catheter
CTP-XXXIP	NuPort-CT, port system incl. introducer set, standard, CT compatible, power injection capable, with PUR catheter
CTP-XXXCS	NuPort-CT, port system set, standard, CT compatible, power injection capable, with silicone catheter
CTP-XXXIS	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter
CTP-L-XXXCP	NuPort-CT, port system set, CT compatible, power injection capable, with PUR catheter, left chest placement
CTP-L-XXXIP	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with PUR catheter, left chest placement
CTP-L-XXXCS	NuPort-CT, port system set, CT compatible, power injection capable, with silicone catheter, left chest placement
CTP-L-XXXIS	NuPort-CT, port system incl. introducer set, CT compatible, power injection capable, with silicone catheter, left chest placement
LPA-XXXCP	NuPort-LP, port system set, CT compatible, power injection capable, with PUR catheter
LPA-XXXIP	NuPort-LP, port system incl. introducer set, CT compatible, power injection capable, with PUR catheter
LPA-XXXCS	NuPort-LP, port system set, CT compatible, power injection capable, with silicone catheter
LPA-XXXIS	NuPort-LP, port system incl. introducer set, CT-compatible, power injection capable, with silicone catheter
SPB-XXXCP	C-Port, port system set, standard, with PUR catheter
SPB-XXXIP	C-Port, port system incl. introducer set, standard, with PUR catheter
SPB-XXXCS	C-Port, port system set, standard, with silicone catheter
SPB-XXXIS	C-Port, port system incl. introducer set, standard, with silicone catheter
PH-XXIS-YY	NuPort-CT, port system incl. introducer set, standard, with silicone catheter
LPA-XXXIP-YY	NuPort-LPA, port system incl. introducer set, with PUR catheter
CTP-XXXCP-YYY	NuPort-CT, port system set, CT compatible, with PUR catheter
CTP-XXXIP-YYY	NuPort-CT, port system incl. introducer set, CT compatible, with PUR catheter
CTP-L-XXXCP-YYY	NuPort-CT, port system set, CT compatible, with PUR catheter, left chest placement
CTP-L-XXXIP-YYY	NuPort-CT, port system incl. introducer set, CT compatible, with PUR catheter, left chest placement

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8086551	Fresenius C-Port HF 8F, 8F silicone catheter, with introducer set
8086561	Fresenius C-Port HF 9.6F, 9.6F silicone catheter, with introducer set

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- Esame e valutazione di osservazioni, rapporti, eventi e incidenti relativi ai prodotti PHS Medical GmbH