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# Port catheter care

## index

Who is allowed to puncture port catheters?	6
What is a NuPort/C-Port?	6
How does the handling of port catheters of the type NuPort/C-Port differ from port catheters of other manufacturers?	6
Does the NuPort/C-Port have a valve function?	6
What is the total priming volume of a NuPort/C-Port?	6
Why must syringes with a nominal volume of less than 10 ml not be used for injections via the port catheter?	6
What are the different types of port cannulas?	6
Are there any restrictions on the use of certain port cannulae?	6
Which cannula diameter should be used for which applications?	6
What cannula length should be used?	6
In which way is the NuPort/C-Port to be punctured	7
What may be the cause when patients report pain during an infusion?	7
What to do if the port catheter is not patent or resistance is felt during the injection?	7
When should the port system be flushed?	7
What should be observed after the infusion of drugs?	7
What should be considered during and after infusion of parenteral nutrition solution?	7
What is the procedure for an infusion of blood products?	7
Wie ist das Vorgehen bei und nach einer Blutentnahme?	7
How should the NuPort/C-Port be flushed?	8
How long may port cannulas be used uninterrupted („dwell time“)?	8
What to do in case of infusion and aspiration difficulties?	8
How should the puncture needle be removed from the port?	8
What should be done if air accidentally gets into the port catheter?	8
What to do if you still have questions regardless of these port catheter care and use instructions?	8

PHS Medical GmbH  
Ederweg 3  
34277 Fuldabrück . Germany  
Tel. +49 561 998597-0 . Fax -199  
info@phs-medical.de  
www.phs-medical.com

# Port catheter care and usage information

**Note:** Your institution's recommendations and regulations for the care and use of the port catheter may differ from the recommendations and instructions in this document. In this case, the responsible physician at your facility must decide which procedures to use.

---

## Who is allowed to puncture port catheters?

Handling port catheters is a medical activity which – like other medical activities – can be delegated to trained, medical personnel.

---

## What is a NuPort/C-Port?

The NuPort/C-Port product families are port catheters from the manufacturer PHS Medical GmbH. With their flow-dynamically optimized port chambers, compared to port catheters of other manufacturers, the procedures to be performed through them are improved by a higher efficiency in flow. This also results in higher flushing efficiency, especially after higher viscosity infusions. The port catheters are also available in variants suitable and approved for high-pressure injections, recognizable by the "CT" suffix to the product name or the heart shape of the port housing.

---

## How does the handling of port catheters of the type NuPort/C-Port differ from port catheters of other manufacturers?

The handling of port catheters of the NuPort/C-Port type is not fundamentally different from that of port catheters of other manufacturers.

---

## Does the NuPort/C-Port have a valve function?

No. All port catheters of the NuPort/C-Port families are without valves and the catheters are distally open. This ensures unobstructed flow and also backflow. Blood sampling is therefore technically possible with our systems.

---

## What is the total priming volume of a NuPort/C-Port?

Typically, a total priming volume of 1.0 ml can be assumed for our port catheters with prethoracic implantation.

---

## Why must syringes with a nominal volume of less than 10 ml not be used for injections via the port catheter?

With syringes smaller than 10 ml nominal volume, there is a risk of producing unacceptable high pressure. Due to the smaller syringe plunger diameters, the pressure generated via the syringe increases while the (hand) force remains constant. Even with 5 ml syringes, it is possible with manual force to pressurize port catheters to their bursting point if the port catheter is not patent.

---

## What are the different types of port cannulas?

Non-coring cannulae must be used with port catheters. In addition to the classic port cannula, i.e. Huber ("Huber needle"), there are also the spoon cut form derived from the Huber cannula from various manufacturers, Trocar cannulae and non-coring Pencil-Point cannulae. In case of doubt, please contact the representative of the respective cannula manufacturer responsible for you.

---

## Are there any restrictions on the use of certain port cannulae?

For standard port catheters, all port cannulae on the market with outer diameters up to 19G (1.1 mm) can be used.

**CAUTION:** Port cannulae with diameters larger than 20G must not be used for our low-profile port. Due to the significantly smaller septum area, compared to standard ports, there is a risk of septum damage with larger cannula diameters. For reasons of procurement optimization, we have

observed in the market that for applications via port catheters, only 19G port cannulae are often used in thrapy facilities. When using 19G cannulae in connection with small septum surfaces, there is a risk of septum damage if the port cannulae are not optimally secured due to the longer cutting edge at the cannula outlet. In order to ensure a high level of application and patient safety, this has prompted us to limit the maximum permitted outer cannula diameter for our low-profile ports to 20G (0.9 mm). This significantly reduces potential damage to the port septum. Please contact us in case of doubt and for further questions.

---

## Which cannula diameter should be used for which applications?

The usual cannula diameters for use with port catheters are as follows:

- 22G (ca. 0,7 mm)
- 20G (ca. 0,9 mm)
- 19G (ca. 1,1 mm)

The following is recommended for the selected cannula outside diameters:

- 22G for aqueous infusates,
- 20G for medium viscosity infusates and for blood sampling via low-profile port catheters,
- 19G for blood collection from standard port catheters to further reduce possible hemolysis, for administration of blood products and for infusion of high viscosity parenteral nutrition solution.

**CAUTION:** In this context, please also note our information under "Are there any restrictions on the use of certain port cannulae?".

---

## What cannula length should be used?

The appropriate cannula length must be determined individually for each patient.

**The following should be considered:** The puncture procedure is completed when the puncture needle touches the port floor

- Make sure that the patient's skin is not compressed in the area of the puncture site. In this case, there is a risk that the cannula will change position over the course of the application and the outlet of the cannula will migrate back into the septum. This can cause decreased flow and even complete occlusion of the cannula.

If possible, use a different cannula length in this case or check the position of the puncture cannula more frequently.

- However, the base plate of the port cannula should not be more than approx. 5 mm from the skin surface, as otherwise the fixation of the port cannula is made more difficult. Over the course of the application, this can lead to cannula dislocation. Likewise, the use of cannulae that are too long leads to higher shear forces in the septum, which can damage the septum.

If possible, use a different cannula length in this case. If necessary, you can bridge and stabilize an excessively large gap between the skin surface and the port cannula with an intermediate layer of a sterile gauze compress that has been cut into beforehand. In this case, check the position of the puncture cannula more frequently.

It is recommended to document at least the used cannula length, better also its outer diameter, and type in the patient records.

>>

---

## In which way is the NuPort/C-Port to be punctured

Please proceed as follows for puncturing:

- First carry out the skin disinfection. Observe the exposure and evaporation time of the skin disinfectant.
- Localize the port housing under the skin surface by palpation.
- Grasp the port between the thumb and index finger of one hand and hold it in place.
- Now puncture the port/port septum with the other hand using a port cannula suitable for this purpose and previously rinsed with saline.
- Insert the cannula vertically through the skin and septum without tilting or circular motion using gentle pressure until the cannula touches the port chamber floor. Here, it is important to press the cannula tip only with very light pressure on the port chamber floor, as the tip can be deformed if too much pressure is applied. This may result in permanent damage to the septum when the cannula is withdrawn from the port.
- Check the correct cannula position in the port chamber by slow aspiration.
- To confirm correct cannula position and to confirm aspiration function, blood should become visible in the syringe after aspiration of approximately 1 to 1.5 ml without resistance.
- Discard this aspirate obtained and immediately flush with at least 2 x 10 ml saline solution.

Now you can commence infusion of the intended drug or infusate.

**CAUTION:** Reuse of a port cannula (e.g., on the same patient for puncture correction) is not recommended due to the risk that a cannula tip bent or torn during the initial puncture will permanently damage the septum if used again.

**CAUTION:** Appropriate aseptic techniques and sterile syringes must be used for each application.

---

## What may be the cause when patients report pain during an infusion?

Possible causes of pain during an infusion can include, for example

- an incorrectly placed port cannula (malpuncture),
- a damaged septum resulting in extravasation,
- a disconnected catheter at the catheter connection of the port housing or
- a defective catheter (catheter rupture).

**CAUTION:** If malpuncture can be ruled out, discontinue use and inform the responsible physician.

---

## What to do if the port catheter is not patent or resistance is felt during the injection?

If you experience resistance during an injection or infusion, stop the procedure immediately. Excessive injection pressure can cause the catheter or port to burst.

Resistance can be caused by

- an incorrectly placed port cannula (malpuncture),
- a dislocated or obstructed catheter lumen, or
- a port occlusion.

Please also see section "What to do in case of infusion and aspiration difficulties?"

**CAUTION:** If malpuncture can be ruled out, discontinue use and inform the responsible physician. Restoration of patency is a physician's activity.

---

## When should the port system be flushed?

The port system must be flushed after each use, for example

- after each injection of medication,
- after the injection of parenteral nutrition solution,
- after blood sampling and
- and, if necessary, should be flushed during longer breaks between therapies.

---

## What should be observed after the infusion of drugs?

After each drug injection, flush the port system with 10 ml of normal saline solution.

---

## What should be considered during and after infusion of parenteral nutrition solution?

Due to the wide range of composition of parenteral nutrition solutions, there are different viscosities which can lead to different and possibly longer infusion times. (see also "Which cannula diameter should be used for which applications?"). It is important to flush the port catheter thoroughly with saline solution immediately after completing an infusion of parenteral nutrition solution. Experts recommend flushing using the push-pause technique with at least 2 x 10 ml of saline.

---

## What is the procedure for an infusion of blood products?

After infusing blood product, immediately flush the port catheter with at least 2 x 10 ml saline solution using the push-pause technique. Flushing with 3 x 10 ml can improve the flushing result.

**CAUTION:** It is extremely important to flush the port with sufficient saline solution immediately after infusing blood products. If blood products remain in the port catheter for a longer period of time, they may clot and lead to occlusion of the port catheter.

---

## What is the procedure during and after a blood collection?

For blood collection, the following procedure as well as the associated precautions must be observed:

1. Insert the port cannula (see also: „Which cannula diameter should be used for which applications?“) at the prepared site. Check the patency of the system by slow aspiration, and then flush with 20 ml of saline solution.
2. Take a sample of at least 5 ml of blood and discard it together with the syringe. Use a new port cannula for a new puncture. (see also: „What to do in case of infusion and aspiration difficulties?“).
3. Connect syringes or collection tubes as needed to obtain the necessary blood samples. Aspirate using slow aspiration to avoid hemolysis due to negative pressure.
4. Immediately flush the port catheter after the last blood collection, using the push-pause technique, with at least 2 x 10 ml saline.

**CAUTION:** It is extremely important to flush the port with sufficient saline solution immediately after blood collection. If blood remains in the port catheter for a longer period of time, this may clot and lead to occlusion of the port catheter.

>>

### How often should the NuPort/C-Port be flushed during therapy intervals?

PHS Medical recommends flushing the ports of the NuPort/C-Port families with saline solution every three months during longer intervals between applications to prevent occlusion and to ensure that the port is still functional. However, the frequency of flushing is a decision based on the type of use. Here, flushing intervals of once a week, once a month, every three months, or even longer intervals are common in practice. Taking into account that every port puncture is also associated with the risk of germ introduction, the intervals should be determined individually.

### How should the NuPort/C-Port be flushed?

It has been shown that the use of the so-called "push-pause" technique improves the flushing effect. A push-pause rhythm with intervals of approx. ½ second, applied over the entire duration of the irrigation, is advantageous for the flushing effect. Due to the induced pressure pulses, adhesions to the inner walls of port catheters are better dislodged and thus the overall flushing effect is improved.

### How long may port cannulas be used uninterrupted ("dwell time")?

There is no uniform expert recommendation on the duration of use of port cannulas. Please follow the standard applicable in your institution for the maximum dwell time of port cannulas.

### What to do in case of infusion and aspiration difficulties?

Infusion and aspiration difficulties can have a variety of causes such as.

- incorrect puncture of the port,
- unfavorable patient positioning,
- patient is agitated resulting in high muscle and vascular tone,
- catheter kinking or entrapment (pinch-off syndrome),
- distal catheter dislocation,
- distal sidewall position of the catheter,
- fibrin exudation of the catheter with formation of a pseudovalue,
- biofilm exudation of the intravascular catheter segment, or
- port occlusion

Improved patency or flow may be achieved by postural modification involving the shoulder girdle.

**CAUTION:** It is advisable to inform the responsible physician about the infusion or aspiration difficulties and to document the situation.

**CAUTION:** If malpuncture can be ruled out, discontinue use and inform the responsible physician. Restoration of patency is a physician's activity.

### How should the puncture needle be removed from the port?

When removing the cannula from the port, the flushing process should be continued at a lower flow rate. This is done by further reducing the pressure on the syringe plunger when removing the cannula. This means that the flushing process is continued with reduced flow during the removal of the cannula.

**CAUTION:** When the port cannula is removed from the septum of the port, a lifting movement occurs in the septum. This lifting movement creates a negative pressure, which aspirates blood at the catheter tip. This can later lead to occlusion of the catheter. This is safely prevented by generating a flow when removing the cannula.

### What should be done if air accidentally gets into the port catheter?

**CAUTION:** Please contact the responsible physician immediately.

### What to do if you still have questions regardless of these port catheter care and use instructions?

You are welcome to contact us by telephone or e-mail with any questions you still have. In this way, we can hopefully help you in the short term and also support other users by supplementing these notes. We look forward to a contact and dialog with you.

PHS Medical GmbH  
Ederweg 3  
34277 Fuldabrück . Germany  
Tel. +49 561 998597-0 . Fax -199  
info@phs-medical.de  
www.phs-medical.com